
Instructions for Use

Prep-Lock™

Tamper Evident Cap for Luer Lock I.V. Syringes

en	Instructions For Use
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
no	Bruksanvisning
sv	Bruksanvisning
es	Instrucciones de uso
fr	Instructions d'utilisation
it	Istruzioni per l'uso

Tamper Evident Cap for Luer Lock I.V. Syringes	en
Originalitätsverschlusskappe für i.v. Luer-Lock- Spritzen	de
Verzegelde dop voor Luer-lock IV. Spuiten	nl
Berøringssikret hette for intravenøse Luer-lock-sprøyter	no
Plomberat lock för intravenös luerlock. Sprutor	sv
Tapa con sensor de apertura para Luer Lock para administración por vía i.v. Jeringuillas	es
Capuchon inviolable pour seringues intraveineuses avec raccord Luer Lock	fr
Tappo antimanomissione per siringhe Luer Lock IV	it

Emergo Europe
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Emergo Australia
201 Sussex Street
Darling Park, Tower II, Level 20
Sydney NSW 2000

Emergo UK
Compass House Vision Park Chivers Way
Histon Cambridge, CB24 9BZ
United Kingdom

MADE IN THE USA
International Medical Industries, Inc.
2981 Gateway Dr.
Pompano Beach, FL 33069 USA
(954) 917-9570
U.S. Patent: imiweb.com/patents

Instructions For Use

Tamper Evident Cap for Luer Lock I.V. Syringes

Indications for use

- The Tamper Evident Cap is a single-use device intended to provide evidence of use, positive closure and maintain sterility for Luer lock Syringes.

Operating Instructions

Use Aseptic Technique through all steps.

- Hold package firmly at tab and peel open while maintaining aseptic technique.
- Insert IV syringe's male luer into the cap and rotate clockwise while pushing downwards.

Verify the installation by rotating the cap until snug.

- Remove the IV syringe and Tamper Evident Cap from package.
- To remove the Tamper Evident Cap, grip the Tamper Evident Cap and pull the outer sleeve straight off.
- Twist off the tip cap counter clockwise and remove the small ring.

Damaged components or missing cap indicates potential tampering.

Contraindications

- None when used as described in the operating instructions.

Cautions / Warnings

- Do not use if sterile barrier is broken.

Potential adverse effects if device is re-used

- Cross contamination

Intended users

- Healthcare practitioners

Patient groups

- General hospital

Disposal

- After use, discard the device following the standard procedures for hospital waste.

Materiovigilance and Traceability

- This product is within the scope of materiovigilance. Any incident or side effect, associated, or suspected to be associated with the use of this product must be reported to the manufacturer and competent authority as listed on the product label.

Symbols Applying To This Device

	Unique Device Identifier
	For sale by or on the order of a licensed healthcare practitioner
	Date of Manufacturing / Country of Manufacture
	Keep dry
	Keep away from sunlight
	Sterilized using ethylene oxide
	Do not re-use

	Use-by-date
	Caution
	Do not re-sterilize
	Do not use if package is damaged
	Single sterile barrier
	Non-pyrogenic
	Consult instruction for use

	Medical Device
	Catalogue number
	Batch code
	Authorized Representative in the European Union
	Manufacturer

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