

# VETERINARY CARE PRODUCTS





**04** Medication Delivery

**24** Bags

**28** Boxes & bins

**30** Consumer Products

**34** Security Seals

**35** PPE

**40** Labels



## PLASTIC MEDICINE BOTTLES, AMBER

### Amber Plastic Oval Medicine Bottles with Caps

- USP Compliant: Child-resistant cap conforms to USP "Tight" standards.
- Better Accuracy: Easy-to-read graduations
- Label Ready: Flat back panel for easy label placement
- UV Protection: Ideal for light-sensitive liquid medication

ITEM NO.	SIZE	BOTTLE OPENING	CASE QTY
19489	2 oz.	24 mm	100
19490	4 oz.	24 mm	100
19491	8 oz.	24 mm	100
20833	16 oz.	28 mm	50



## GLASS BOTTLES, AMBER

### Amber Glass Bottles & Accessories

- UV Light Safe - Amber Type 3 pharmaceutical glass bottle is ideal for dispensing topical treatments and other UV light-sensitive liquid medications.
- Flexible - Bottles provide usage flexibility with optional accessories such as glass droppers (19456-19461) and tamper evident caps (10067, 10068).
- Multiple Sizes - Choose from six container sizes to fit your pharmaceutical packaging and dispensing needs.

ITEM NO.	SIZE	CASE QTY	DROPPER CAP
19450	5 ml	340	19456
19451	10 ml	255	19457
19452	20 ml	216	19458
19453	30 ml	168	19459
19454	50 ml	99	19460
19455	100 ml	68	19461



## CUPS

### Plastic Cups

- Available in a variety of sizes
- Available graduated

ITEM NO.	SIZE	GRADUATED	CASE QTY
7364	1 oz.	YES	1000
7362	2 oz.	YES	1000
7360	3 oz.	NO	100
7361	5 oz.	NO	50




## DISPENSING PIN

### ExactaMed® Dispensing Pin

- Ideal for vented additions and/or withdrawals of fluid from multi-dose vials
- Molded plastic spike with filtered vent, sheath and dust cap
- Eliminate risk of air being drawn into solution - prevents air bubbles from being aspirated back into the syringe; with offset air inlet, it lets air in above the fluid
- Pin housing and screw caps are made out of acrylonitrile butadiene styrene; Sheath is made out of polyvinyl chloride; vent assembly housing is made out of polypropylene with a 3 micron acrylic copolymer membrane - hydrophobic and oleophobic vent
- Sterile and non-pyrogenic in unopened, undamaged package
- Latex free
- 50 per package, packaged individually

ITEM NO.	CASE QTY
47	50



## DROPPER BOTTLES, STERILE

### Plastic Dropper Bottles, Sterile

- Sterility Assurance - Ideal for the preparation of sterile ophthalmic solutions. Bottle exterior is sterile until the package is opened; the inside is sterile until the cap is removed.
- All-in-One System - Includes snap-in dropper tips and screw on caps. Releases approximately .028mL per drop.
- Radiation sterilized
- USP 797 Compliant

ITEM NO.	SIZE	AMOUNT RELEASED PER DROP	CASE QTY
7780	3 mL	0.0238 ml/g	12
7780-02	3 mL	0.0238 ml/g	1
7782	7 mL	0.035 ml/g	12
7782-01	7 mL	0.035 ml/g	1
7784	15 mL	0.035 ml/g	12
7784-01	15 mL	0.035 ml/g	1
7786	30 mL	0.035 ml/g	12
7786-01	30 mL	0.035 ml/g	1



## PLASTIC BOTTLES & ACCESSORIES

### Cylinder Plastic Bottles

- Versatile - Bottles can be used for storing, dispensing, mixing, sampling, etc.
- Cap Options - Child-resistant screw caps, flip-top caps and spray pumps are sold separately. See options below
- HDPE bottles have 20.6mm openings

ITEM NO.	SIZE	CASE QTY
10407	100 mL	50
10408	250 mL	54
10409	500 mL	28
10410	1000 mL	60



### Accessories

- Flip-Top Caps: Great for pouring without spilling and dispensing product by squeezing the bottle
- Child-Resistant Screw Caps: Provide tamper-evident closure for Cylinder Plastic Bottles
- Versatile: One size cap fits all bottle sizes

ITEM NO.	DESCRIPTION	CASE QTY
10411	Screw Caps	100
10412	Child Resistant Screw Caps	100
10413	Flip top Cap	100



## ORAL DISPENSERS, CLEAR

### ExactaMed Oral Dispensers, Clear

- Self-righting tip caps included (will not fit luer needle)
- Safe for neonatal, pediatric and adult use
- Unique tip design prevents wrong-route administration
- Chemically inert to provide stability for drug storage and transport
- Blue printing, clear labeling and white plunger provide distinct differentiation from IV syringes
- Easy to read dose markings, even when filled with dark liquids
- Enlarged numbers for quick and accurate volume identification
- Includes distinct bilingual labeling “For Oral Use Only”
- Latex-Free

ITEM NO.	SIZE	INCREMENTS	QTY
10200	0.5 ml	0.01 ml	100
200	0.5 ml	0.01 ml	500
10201	1 ml	0.01 ml	100
201	1 ml	0.01 ml	500
10203	3 ml	0.10 ml	100
801	1 mL	0.01 mL	2000
901	1 mL	0.01 mL	3000
203	3 ml	0.01 ml	500
803	3 ml	0.01 ml	2000
10205	5 ml	0.2 ml	100
205	5 ml	0.2 ml	500
605	5 ml	0.2 ml	1500
805	5 mL	0.2 mL	2000
10210	10 mL	0.2 mL	100
210	10 mL	0.2 mL	500
710	10 mL	0.2 mL	1800
810	10 mL	0.2 mL	2000
10220	20 mL	1 mL	100
1221	60 mL	1 mL	75



## ORAL DISPENSERS, AMBER

### ExactaMed Oral Dispensers, Amber

- Amber dispensers are ideal for dispensing light-sensitive drug applications
- Self-righting tip caps included (will not fit luer needle)
- Safe for neonatal, pediatric and adult use
- Unique tip design prevents wrong-route administration
- Chemically inert to provide stability for drug storage and transport
- White printing, clear labeling and white plunger provide distinct differentiation from IV syringes
- Easy to read dose markings, even when filled with dark liquids
- Enlarged numbers for quick and accurate volume identification
- Includes distinct bilingual labeling “For Oral Use Only”
- Latex-Free

ITEM NO.	SIZE	INCREMENTS	QTY
11200	0.5 mL	0.01 mL	100
1200	0.5 mL	0.01 mL	500
11201	1 mL	0.01 mL	100
1201	1 mL	0.01 mL	500
1801	1 mL	0.01 mL	2000
1901	1 mL	0.01 mL	3000
11203	3 mL	0.01 mL	100
1203	3 mL	0.01 mL	500
1803	3 mL	0.01 mL	2000
1903	3 mL	0.01 mL	2500
11205	5 mL	0.2 mL	100
1205	5 mL	0.2 mL	500
11605	5 mL	0.2 mL	1500
1805	5 mL	0.2 mL	2000
11210	10 mL	0.2 mL	100
1210	10 mL	0.2 mL	500
11710	10 mL	0.2 mL	1800
1810	10 mL	0.2 mL	2000
11220	20 mL	1 mL	100
1222	60 mL	1 mL	25



## ADAPTA CAPS

- Easy to Use: Converts a standard medication bottle into a filling device for oral dispensers
- Material: Low-density polyethylene
- Available in multiple sizes
- Quantity: 25 per package or 100 per case

### Adapta Caps

ITEM NO.	BOTTLE NECK SIZE	PACK QTY	CASE QTY
51001	A-18 mm	25	100
51002	B-20 mm	25	100
51004	D-23 mm For glass bottles	25	100
51005	E-28 mm Short neck	25	100
51006	F-28 mm Long neck	25	100
51013	M-24 mm For plastic bottles	25	100

MOST POPULAR SIZE!



## PRESS-IN BOTTLE ADAPTERS

- Mess-Free Filling - Use Press-In Bottle Adapters with liquid medication bottles to easily fill oral dispensers without creating a mess
- Hassle-Free Design - Prescription bottle caps can be used with plugs in place. Adapters provide a no-leak closure
- Size Options - Please note that the bottle adapter insert sizes listed correspond to the outer neck of bottle and not the inside opening

ITEM NO.	FOR BOTTLE NECK	QTY
17282	20 mm	20
17282-10	20 mm	100
17283	24 mm	20
17283-10	24 mm	100
17284	28 mm	20
17284-10	28 mm	100
17285	33 mm	20
17285-10	33 mm	100
19392	22mm	20
19392-10	22mm	100
19393	24mm	20
19393-10	24mm	100
19394	28mm	20
19394-10	28mm	100

MOST POPULAR SIZES!



## STERILE RAPID FILL CONNECTORS

### Rapid Fill Connector, Luer Lock to Luer Lock with Cap

- For filling sterile dispensers from a luer connection
- Luer lock-to-oral slip connection
- Connects to male luer or syringes, tube sets, repeating syringes, bag ports, etc.
- Flange facilitates contamination-free handling
- Sterile
- Made out of natural polypropylene
- 50 per package, packaged individually

ITEM NO.	DESCRIPTION	COLOUR	QTY
12703	Luer Lock-to-Oral Slip	Natural	50
13701	Luer Lock-to-Luer Slip	Blue	50
13901	Luer-Lock-to-Luer Lock	Red	50
13915	Luer Lock-to-Luer Lock with Cap	Red	50



## TAMPER TAPE

### Integrity Tamper-Evident Tape Strips | 1281

- Compatible: Tape strip molds easily to almost any item.
- Tape reads: Integrity Seal, Sealed by and Date.
- Quantity - 100 strips per box
- Dimensions - 6 $\frac{5}{8}$ " L x 1 $\frac{3}{8}$ " W

### Tamper Resistant Tape

- Convenient: Perfect for use over items such as jar tops, bottle caps, syringes, boxes, specimen vials and similar containers that are hard to secure with ordinary security seals to help deter and show tampering. Tape molds easily to almost any item.
- Deters Tampering - Transparent Red polystyrene tape with serrated edges is printed with the message: "TAMPER RESISTANT TAPE" in Black and will tear in small pieces to show evidence of tamper upon removal.
- Easy to Use - Wax paper liner is slit lengthwise and allows for easy removal from roll.

ITEM NO.	DESCRIPTION	DIMENSIONS	COLOUR	QTY
12281	Integrity Tamper-Evident Tape Strips	6 $\frac{5}{8}$ " L x 1 $\frac{3}{8}$ " W	Red	100 Strips
8000	Tamper Resistant Tape	108'L x 1/2"H	Red	2 Rolls
8002	Tamper Resistant Tape	108'L x 1/2"H	Clear	2 Rolls
8008	Tamper Tape	108'L x 1/2"H	Orange	2 Rolls



## TAPE & ACCESSORIES

### Handy Padlock Seals, Unnumbered

- Quick Access - Handy padlock seal allows quick access to items in crash cart boxes, on carts and in other spaces where products need security.
- Average breaking strength - 10-13 lbs.
- Material - HDPE and Polypropylene
- Tool needed for removal - No
- Autoclavable - Up to 325°F (163°C)
- Colors - Black, Blue, Green, Red, White and Yellow. Specify color when ordering.
- Quantity - 100 per package
- Dimensions - 1-3/8"W x 1-5/8"H



### Tamper-Evident Syringe Seals, Blank

- Multi-Use: Non-sterile seals can be used on all syringes, vials, boxes, otic/ophthalmic containers, etc.
- Molds Easily : Flexible foil seals attach to syringe caps or sheaths. Cannot be reapplied once removed.
- Colors: Blue, Red and Yellow. Specify color when ordering.

ITEM NO.	TAPE WIDTH	TAPE LENGTH	CORE SIZE	QTY
T1-1	1" (24mm)	33 yards	1"	1 Roll
T1-3	1" (24mm)	72 yards	3"	1 Roll
T1-5	1.5" (38mm)	72 yards	3"	1 Roll
T1-7	2" (50mm)	55 yards	3"	1 Roll
TD1	Tape Dispenser For 1" Wide Tape (1" & 3" Cores)			1
GF2306	Tape Dispenser For 2" Wide Tape			1
GF2336	Tape Dispenser For 3" Wide Tape (3" Cores)			1

## BEAKERS

### Beakers

- Scale: White enamel
- Clean Results: Spout for easy pouring
- Compatible: Safe for use on hot plates

ITEM NO.	VOLUME	INCREMENTS	QTY
3372	30 mL	10 mL	1
3024	50 mL	10 mL	1
3025	100 mL	10 mL	1
3373	150 mL	25 mL	1
3026	250 mL	25 mL	1
3374	400 mL	25 mL	1
3027	600 mL	50 mL	1
3375	1000 mL	50 mL	1



## CYLINDERS

### Glass Graduated Cylinder

- Made from high-quality borosilicate glass
- Single metric scale
- Has an easy-pour spout with large circular base
- Rated at 20°C
- Not recommended for autoclaving; for autoclaving please refer to 1634AA series

ITEM NO.	VOLUME	INCREMENTS	QTY
3002	10 mL	0.2 mL	1
3009	25 mL	0.5 mL	1
3006	50 mL	1 mL	1
3008	100 mL	1 mL	1
3014	250 mL	2 mL	1
3015	500 mL	5 mL	1
3371	1000 mL	10 mL	1



## PILL COUNTERS


### Fretwell's Triangle

- Precision-made scientific instrument
- Count regular shaped tablets with absolute accuracy and speed
- Simply count the number of rows to determine tablet count - unique legend determines count
- Tablets glide easily on the stainless steel triangle
- Easy to clean

### Pill Counters

- Plastic pill counting tray

ITEM NO.	DESCRIPTION	QTY
FT1	Fretwell's Triangle	1
D-50	Pill Counter With Spatula	1
7147	Tablet Counter, Right-Handed	1
7148	Tablet Counter, Left-Handed	1



## MEDICATION ENVELOPE

### Glass Graduated Cylinder

- Cost Effective - Sealable envelope for small quantities of medication
- Write-On Area - For patient information and instructions

ITEM NO.	DESCRIPTION	DIMENSIONS	QTY
8020	Medication Envelopes	3½"W x 2¼"H	250



## MORTAR & PESTLES

### Flint Glass Mortar & Pestles

- Convenient - Use for preparing medications for patients or compounding recipes
- Easy to Clean - Can be sterilized, cleaned with alcohol or washed with soap and water
- Note: Not recommended for use with drastic temperature changes
- Flint glass models are ideal for compounding liquids, ointments & pastes

ITEM NO.	DESCRIPTION	SIZE	QTY
3070	Flint Glass Mortar	2 oz. (60 mL)	1
3080	Flint Glass Pestle	2 oz. (60 mL)	1
3071	Flint Glass Mortar & Pestle Set	2 oz. (60 mL)	1
3072	Flint Glass Mortar	4 oz. (120 mL)	1
3082	Flint Glass Pestle	4 oz. (120 mL)	1
3073	Flint Glass Mortar & Pestle Set	4 oz. (120 mL)	1
3074	Flint Glass Mortar	8 oz. (240 mL)	1
3084	Flint Glass Pestle	8 oz. (240 mL)	1
3075	Flint Glass Mortar & Pestle Set	8 oz. (240 mL)	1
3076	Flint Glass Mortar	16 oz. (480 mL)	1
3086	Flint Glass Pestle	16 oz. (480 mL)	1
3382	Flint Glass Mortar & Pestle Set	16 oz. (480 mL)	1



## COMPOUNDING ACCESSORIES

### Ointment Pad

- 12" square
- 100 sheets per pad

### Stainless Steel Spatula

- Versatile - Flexible stainless steel blade with comfortable wooden handle
- Easy Cleaning - Use soap and water

### Plastic Weighing Boats

- Inexpensive containers for weighing liquids or solids.
- Resilient polystyrene construction
- Non-stick surface, rugged enough for mixing
- Stack and nest easily

ITEM NO.	DESCRIPTION	DIMENSIONS	QTY
3088	Ointment Pad	12" x 12"	100/pad
3042	Spatula Stainless Steel	4" Blade	1
3043	Spatula Stainless Steel	5" Blade	1
3028	Spatula Stainless Steel	6" Blade	1
3029	Spatula Stainless Steel	8" Blade	1
3429	Spatula Stainless Steel	10" Blade	1
3430	Spatula Stainless Steel	12" Blade	1
3420	Plastic Weighing Boat Micro	4.4 x 4.4	100
40100	Plastic weighing Boat, Small	4.5 x 4.5	100
3090	Plastic weighing Boat, Medium	8.3 x 8.3	100
3092	Plastic weighing Boat, Large	14 x 14	100



## VIALS, STERILE

- Sterile Processed: Type 1, USP, borosilicate glass vials are dry-heat sterilized and pyrogen free.
- Compatibility: Empty Sterile Vial stoppers are latex free.
- Please note: Caps are sold separately


### Amber Vials

ITEM NO.	CAPACITY	CAP SIZE	QTY
19812	2 mL	13 mm	100
19813	5 mL	20 mm	25
18702	10 mL	20 mm	25
18703	30 mL	20 mm	25
18704	50 mL	20 mm	25

### Clear Vials

ITEM NO.	CAPACITY	PACK QTY	CASE QTY
7760-10	2 mL	-	100
7761-03	5 mL	25	100
7762-03	10 mL	25	100
7763-10	30 mL	25	100
7765-03	50 mL	25	-
7764-06	100 mL	25	-

→ MOST POPULAR SIZE!



## UV ZIP IT BAGS, AMBER

- Easy to access and provides sturdy storage
- Latex-free
- Protects light sensitive medications - amber and green colour bags absorb 92% to 98% of light in the UV range - the dark brown and black colour bags absorb 99.9% of light in the UV range
- Lighter colours allow the nurse to read a label or check the level of an IV fluid through the UVLI bag
- Slit-top bags provide a professional look



ITEM NO.	SIZE	QTY
AUZ35	3 X 5"	1000
AUZ39	3 X 9"	1000
AUZ46	4 X 6"	1000
AUZ58	5 X 8"	1000
AUZ69	6 X 9"	1000
AUZ814	8 X 14"	1000
AUZ912	9 X 12"	1000

## PRESCRIPTION BAGS & ACCESSORIES

### Hanging Prescription Bags

- Simplify prescription organization and retrieval
- Made of tough, clear 4 mil plastic for extreme durability
- Fits on any rod with a diameter of 1 cm or smaller

### Wall Rack

- Simple solution to prescription organization
- Saves time, improves accuracy and enhances customer satisfaction
- Aluminum frame is 3 ft. wide and holds five 1 cm thick rods
- Each rod holds approximately 8-10 bags per foot or as many as 150 bags per rack
- Easy assembly, all hardware included
- Measures: 36" W x 57.5" H x 7.5" D
- Use with hanging bags and ID tags (sold separately)

ITEM NO.	DESCRIPTION	SIZE	QTY
17540	Hanging Prescription Bag	10.5 x 7.5"	10
17541	Hanging Prescription Bag	12 x 9"	10
17542	Hanging Prescription Bag	14 x 12".5"	10
543	Wall Rack for Prescription Bags	-	1 frame, 5 rods
44	Green ID Tags	-	25



## IV COVERS & BAGS

- UV Blocking - Covers block 95% of UV light and are perfect for dispensing IV admixtures that contain light-sensitive drugs.
- Opacity - IV Covers are Dark Amber
- Transparent - Bags are transparent, but will mask many differences in IV fluid color.
- Thickness - 1.5-mil - 1.8 mil

ITEM NO.	SIZE	DESCRIPTION	TINT	QTY
7589	8 x 14	IV Cover, slit at sealed end	Light	100
7587	6 x 10	IV Cover, slit at sealed end	Amber	100
7587D	6 x 10	IV Cover	Dark	100
7589D	8 x 14	IV Cover, slit at sealed end	Dark	100
7590D	8 x 14	IV Cover	Dark	100
19643	12 x 18	IV Cover	Amber	100
7583	8 x 14	IV Cover, slit at sealed end	Green	100
7584	6 x 10	IV Cover, slit at sealed end	Green	100
7590L	8 x 14	IV Cover	Amber	100
7579	5 x 8.5	IV Cover, slit at sealed end	Green	100



## BAGS, AMBER

### UV Protection Bags, Amber

- Special Features: Bags are resistant to UV light penetration. Ideal for storing or dispensing light sensitive medications.
- Thickness: 3-mil
- Closure: Single-track zipper

ITEM NO.	SIZE	QTY
7567M	2.5 X 9"	100
7567M	2.5 X 9"	1000
7570M	3 X 5"	100
7570M	3 X 5"	1000
7573M	4 X 6"	100
7573M	4 X 6"	1000
7575M	6 X 8"	100
7575M	6 X 8"	1000
7576M	8 X 8"	100
7576M	8 X 8"	1000
7588M	8 X 14"	100
7588M	8 X 14"	1000
7578M	12 X 12"	100
7578M	12 X 12"	1000

→ **MOST POPULAR SIZE!**



## BOXES

### Divider Boxes

- Smart Design - Recessed stacking ledges protect contents of container below
- Durable Plastic - Resistant to most solvents, chemicals, stress cracks
- Autoclavable - 250°F (121°C) for two minutes.
- Optional Accessories - Dividers, lids and label holders are sold separately
- Customizable - Plastic dividers dividable every 1¼"
- Color Options - Blue, Gray, Red or Light Blue. Specify color when ordering.



### Lids

- Lids include lift-holes providing an area to secure the box with security seals (sold separately)
- Semi-clear snap-on lids fit over the top edges of the box
- Clear deluxe lids fit on inside rim of boxes

ITEM NO.	OUTSIDE DIMENSIONS (INCHES L x W x H)	COLOUR
1700	10.875 x 8.25 x 5	● ● ●
1705	16.5 x 10.875 x 6	● ● ●
1710	16.5 x 10.875 x 8	● ● ●
1715	22.375 x 17.375 x 8	● ● ●
1723	16.5 x 10.875 x 3.5	● ● ●
1724	10.875 x 8.25 x 2.5	● ● ●
1725	10.875 x 8.25 x 3.5	● ● ●
1891	16.5 x 10.875 x 4	● ● ●
1892	22.375 x 17.375 x 10	● ● ●
1893	22.375 x 17.375 x 4	● ● ●
1894	22.375 x 17.375 x 6	● ● ●

ITEM NO.	LID TYPE	COMPATABLE WITH
1701	Snap On Lid	1700
1706	Snap On Lid	1891, 1705, 1710
1716	Snap On Lid	1893, 1894, 1715, 1892
1730	Deluxe Lid	1700
1730-01	Deluxe Lid	1724, 1725
1731	Deluxe Lid	1891, 1705, 1710
1731-01	Deluxe Lid	1723
1732	Deluxe Lid	1893, 1894, 1715, 1892

## TILTERS® TILT FRONT BIN

- Modular bins – stack and customize to fit your storage and organization requirements
- Front opens at a 45 degree angle for easy access
- Clear fronts provide easy visual access to inventory
- Moulded in tough beige or white plastic. Specify colour when ordering.

ITEM NO.	OVERALL DIMENSIONS (INCHES L x W x H)	INSIDE DIMENSIONS (INCHES L x W x H)	COMPARTMENTS PER UNIT	COLOUR
1203A	23.625 x 4.5 x 3.625	3.125 x 2.75 x 2.625	6	○ ○
1205E	23.625 x 6.5 x 5.25	4 x 4.5 x 3.75	5	○ ○
1207	23.625 x 8.125 x 6.62	5 x 5.25 x 5	4	○ ○
1210A	23.625 x 9.5 x 7.75	7.75 x 23.625 x 9.5	3	○ ○
1215	23.625 x 13.75 x 9.62	11.875 x 23.625 x 13.875	2	○ ○
1213	11.812 x 13.87 x 12.25	10.5 x 9.375 x 8.25	1	○ ○



## PILL PLANNERS

- Helps plan weekly medication, vitamins and supplements
- Translucent design – easily view contents
- Assorted colours available



ITEM NO.	DESCRIPTION	COMPARTMENTS PER UNIT	QTY
PS050	Mini Pill Boxes	1	2
PS033	Mini Pill Pods	1	3
PS255	Weekly Pill & Planner, Small	7	1
PS265	Weekly Detachable Pill Planner	7	1
PS060	Stackable Pill Planner	7	1
PS280	Pill Planner Twice Daily, XL	14	1
PS118	Weekly Pill Planner With 4x Daily, XL	4	1

## PILL SPLITTERS & CRUSHERS

### Pill Splitter

- Splitter is designed to hold pills securely in place for precise splitting
- Easy to use, safe and convenient
- Blade is made out of surgical stainless steel
- Built-in storage container is designed to catch pills

### Pill Crusher

- Turns tablets into powder effortlessly
- Dual storage compartments
- Cap is threaded for extra security

ITEM NO.	DESCRIPTION	IMPRINT AVAILABLE	QTY
PS7345	Pill & Tablet Splitter	YES	1
PS7325	Deluxe Pill & Tablet Splitter	-	1
PS7346	Pill Crusher	-	1

→ MOST POPULAR SPLITTER

**CUSTOM IMPRINTING AVAILABLE!**  
PERSONALIZE YOUR PRODUCTS TODAY.



## CONSUMER PRODUCTS

### Thermometer | PS720

- For axillary or oral use only
- Large and easy to read liquid crystal display
- Fast and accurate - 10 seconds read time after thermometer reaches a stable temperature reading
- Auto shut-off feature after 10 minutes of inactivity
- Alarm sounds to notify whether temperature is within normal range
- Includes previous temperature memory feature
- Low battery indicator – replaceable battery
- Flexible tip for convenience and safety
- Packaged with probe protection case



### Hot & Cold Compress | HC1400

- Convenient, flexible, reusable gel packs
- Ideal for first-aid application, or keeping items chilled for storage and transport

### Nasal Aspirator | PS900

- Doctor recommended Nasal Aspirator
- Removes nasal congestion
- Use to clean mucus from nasal passageways

ITEM NO.	DESCRIPTION	QTY
PS720	10 Second Flex Tip Digital Thermometer	1
HC1400	Hot & Cold Compress, 7x10"	1
HC4624	Hot & Cold Compress,	
PS900	Nasal Aspirator, 1 Oz	1



## ORAL SYRINGES

### Oral Syringe Imprints

- An ideal way to promote your pharmacy - give-away with filled liquid prescriptions and to new customers
- Include up to 3 lines of text, 30 characters per line, enter details below - contact us for details on how to imprint your company logo
- Imprint in red, blue, green or black ink - custom colours are available by special request (additional pricing may apply)
- Individually poly-bagged
- Graduated with both millilitres and teaspoon measurements
- Latex Free


ITEM NO.	COLOUR	SIZE	IMPRINT AVAILABLE	MIN QTY
IMP205PSW	White Plunger	5 mL	YES	500
IMP205PSB	Blue Plunger	5 mL	YES	500
IMP210PSW	White Plunger	10 mL	YES	500
IMP210PSB	Blue Plunger	10 mL	YES	500



## SEALS


### Handy Padlock Seals, Unnumbered

- Quick Access - Handy padlock seal allows quick access to items in crash cart boxes, on carts and in other spaces where products need security.
- Average breaking strength - 10-13 lbs.
- Material - HDPE and Polypropylene
- Tool needed for removal - No
- Autoclavable - Up to 325°F (163°C)
- Colors - Black, Blue, Green, Red, White and Yellow. Specify color when ordering.
- Quantity - 100 per package
- Dimensions - 1-3/8"W x 1-5/8"H



### Tamper-Evident Syringe Seals, Blank

- Multi-Use: Non-sterile seals can be used on all syringes, vials, boxes, otic/ophthalmic containers, etc.
- Molds Easily : Flexible foil seals attach to syringe caps or sheaths. Cannot be reapplied once removed.
- Colors: Blue, Red and Yellow. Specify color when ordering.



ITEM NO.	DESCRIPTION	DIMENSIONS	QTY
18318	Handy Padlock Seals, Unnumbered	1-3/8"W x 1-5/8"H	100
7868	Tamper-Evident Syringe Seals, Blank	4¼"W x 7/16"H	1000

## SURFACE SAFE™ TWO-STEP APPLICATOR KIT

- Active Ingredients:
- Packet # 1 - High Potency Bleach, contains a towelette in a 2% sodium hypochlorite soap solution.
- Packet # 2 - Bleach Inactivating Agent, contains a towelette in 1% sodium thiosulfate with 0.9% benzyl alcohol.
- On Contact Results - Two-step applicator kit cleans and deactivates chemotherapy drugs on work surfaces to prevent the risk of cross contamination. Creates a cleaner, safer chemotherapy work surface.
- Kit Convenience - Kit consists of two wipes packaged in numbered packets that are applied to the contaminated work surface in numerical order.
- Treatment Area - Each kit treats approximately 2 sq. ft.
- Quantity - 15 kits per box
- Towelette Size - 5 1/2" x 10"

ITEM NO.	DESCRIPTION	SIZE	QTY
18821	Surface Safe™ Two-Step Applicator Kit	5.5" x 10"	15 kits



## BEARD COVERS & BOUFFANT CAPS

### Beard Covers

- Material - Polypropylene
- Sizes - Universal size
- Provides personal protection for doctors, surgeons, nurses and other medical professionals.

### Beard Guards

- Completely cover beards, mustaches, and sideburns.
- Comfortable elastic strap loops behind ears

### Tyvek® Bouffant Caps

- Material - Tyvek material offers superior breathability.
- Construction - Elastic headband provides an excellent barrier for dry particulates including microorganisms and more.

ITEM NO.	DESCRIPTION	QTY
12077	Beard Covers	10
12077-31	Beard Covers	100
HB2060	Beard Guard	100
19346	Tyvek® Bouffant Caps	10



## NITRILE MEDICAL EXAMINATION GLOVES

### Features

- Durable nitrile examination grade glove
- Suitable for various applications such as hospitals, food handling and industrial use
- Textured for extra grip and dexterity
- Excellent tensile strength and touch sensitivity
- Powder free & Latex free
- 4 mil thick
- Blue colour

### Certifications:

- Class 2 MD
- Medical Gloves


ITEM NO.	SIZE	QTY
GNPMS	Small Nitrile Gloves	100/box
GNPMM	Medium Nitrile Gloves	100/box
GNPML	Large Nitrile Gloves	100/box



## FACE MASKS

### 3 Ply Face Mask with Ear Loop

- Elastic Ear Loops, latex free
- Disposable
- Pleated design for maximum comfort
- Lightweight & breathable
- Protection from non-toxic dust and contaminants
- Compliant to EN14683 Standard
- BFE≥98%, PFE ≥98%



### Disposable Full Face Shield

- Construction - Foam headpiece with wide, soft, stretchable strap is lightweight and comfortable. Shield is anti-fogging.
- Dimensions - 7 1/2"L x 13"W



### Features

- Better Comfort - Relieve the tension and irritation caused by elastic bands that rub behind the ears after prolonged periods of face mask use
- Adjustable - Multiple notches on the semi-flexible plastic band allow users to choose their comfort level. Length can also be shortened using scissors.
- Dimensions - 8 5/16"L x 1 7/16"H



ITEM NO.	DESCRIPTION	QTY
B344	3 Ply Face Mask with Ear Loop	50/box
12926	Disposable Full Face Shield	1
20445	Face Mask Ear Relief Bands	10

## GOWNS

### Disposable SMS Isolation Gown

- Disposable gown meets AAMI PB70 Level II barrier protection
- Gown made of non-woven material with full back closure.
- Features: Flat belt, elastic cuffs and tie-on neck closure.
- Colour: Yellow



### Isolation Gowns

- Rated Protection:
  - AAMI PB70:2012 Level 1 and Level 2 protection
  - In compliance with OSHA 1910.1030
  - Passes level 2 AATCC 42 and AATCC 127
- Material - Disposable, latex-free polyethylene gowns offer a cost-effective PPE solution that provides peace of mind for health care workers.
- Construction - Thumb loops offer additional wrist coverage. Slip-on style eliminates tying or fastening around the neck. Open back with tie for comfort. Free of perforations.
- Size - One size fits most



### Disposable Tyvek® Labcoat

- Material - Lightweight Tyvek® material delivers tear strength and puncture-resistance. Low-linting, smooth and opaque. Not water or chemical resistant.
- Construction - Three snaps on the front and elastic wrists.
- Sterilization - Non-Sterile

ITEM NO.	DESCRIPTION	QTY
GW10Y	Disposable SMS Isolation Gown	10
GW100Y	Disposable SMS Isolation Gown	100
20497	Isolation Gowns	50
12084S	Disposable Tyvek® Labcoat	1



# ANESTHESIA/SEDATION LABELS

ITEM NO.	DESCRIPTION	ITEM NO.	DESCRIPTION
B5	Atropine mg/mL	CL-867	Acepromazine
B2	Diazepam mg/mL	CL-706	Acepromazine mg/mL
B31	ePHEDrine mg/mL	CL-997	Alfaxalone
B44	EPIneprine mg/mL	CL-15	Alfaxan
B25	Fentanyl mcg/mL	CL-317	Dexdomitor
B6	Glycopyrrolate mg/mL	CL-1134	Dexmedetomidine
B12	Hydromorphone mg/mL	CL-999	Dexmedetomidine
B41	Ketamine mg/mL	CL-866	Dexmedetomidine mg/mL
B4	Lidocaine mg/mL	CL-1107	Dorminal
B46	Remifentanilmcg/mL	CL-884	Ket/Val...
LPH-87	Bupivacainemg/mL	CL-802	Ketamine
VET40	Butorphanol 10 mg/mL	CL-574	Ketamine
VET27	Dexmedetomidine 500 mcg/mL	CL-33	Ketamine
VET22	Propofol 10 mg/mL	CL-799	Ketamine 50 mg/mL
CL-775	Acepromazine	CL-1114	Romifidine
CL-900	Acepromazine	CL-1106	Sedation

# SURGERY/RECOVERY LABELS

ITEM NO.	DESCRIPTION	ITEM NO.	DESCRIPTION
CL-494	Do Not Flush	CL-731	Medication Added
17573	Heparin Flush	23-2	Medication Added
LPH-21	Hep-Saline	CL-199	Medication Added
CL-708	Heparine/Saline ___ UI/ml	2403	Medication Added to I.V.
MAR1	Medication Added...	LPH-64	Normal Saline
21-2	Medication Added...	CL-1044	Normal Saline
MAR3	Medication Added...	CL-1131	Patient Name...
CL-509	Medication Added		

**Medication Added**  
 Drug additive \_\_\_\_\_  
 Concentration of drug added \_\_\_\_\_  
 Quantity added to bag \_\_\_\_\_  
 Solution concentration in bag \_\_\_\_\_ per ml  
 Added by \_\_\_\_\_  
 Date \_\_\_\_\_ Time \_\_\_\_\_  
 Exp. Date \_\_\_\_\_  
 Patient name \_\_\_\_\_

**MEDICATION ADDED**  
 DRUG \_\_\_\_\_  
 AMOUNT \_\_\_\_\_  
 ADDED BY \_\_\_\_\_  
 DATE \_\_\_\_\_ TIME \_\_\_\_\_  
 EXP. DATE \_\_\_\_\_  
THIS LABEL MUST BE AFFIXED TO ALL INFUSION FLUIDS CONTAINING ADDITIONAL MEDICATION

**MEDICATION - ADDED**  
 PATIENT \_\_\_\_\_ ROOM NO. \_\_\_\_\_  
 DRUG \_\_\_\_\_ DOSAGE \_\_\_\_\_  
 DATE \_\_\_\_\_ TIME \_\_\_\_\_ ADDED BY \_\_\_\_\_  
 EXPIRES DATE \_\_\_\_\_ TIME \_\_\_\_\_

**MEDICATION ADDED TO I.V.**  
 PTS. NAME \_\_\_\_\_  
 ADDED BY/DATE \_\_\_\_\_  
 RETURN TO PHARMACY AFTER \_\_\_\_\_ HOURS

**NORMAL SALINE**

**NORMAL SALINE**

Patient name: \_\_\_\_\_  
 Drug additive: \_\_\_\_\_  
 Concentration of drug added: \_\_\_\_\_  
 Quantity added to bag: \_\_\_\_\_  
 Solution concentration in bag: \_\_\_\_\_ per ml  
 Added by: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 Time: \_\_\_\_\_  
 Exp. Date: \_\_\_\_\_



## GENERIC DRUG LABELS

ITEM NO.	DESCRIPTION	ITEM NO.	DESCRIPTION
VET10	Midazolam 5 mg/mL	VET26	Acepromazine ___mg/mL
VET11	Diazepam 5 mg/mL	VET27	Dexmedetomidine 500 mcg/mL
VET12	Meloxicam 5 mg/mL	VET28	Epinephrine 1:1000
VET13	Atropine 0.5 mg/mL	VET29	Metronidazole 5 mg/mL
VET14	Glycopyrrolate 0.2 mg/mL	VET30	Ampicillin 100 mg/mL
VET15	Fentanyl 50 mcg/mL	VET31	Enrofloxacin 50 mg/mL
VET16	Remifentanyl ___ mg/mL	VET32	Famotidine 10 mg/mL
VET17	Rocuronium 10 mg/mL	VET33	Hydromorphone 10 mg/mL
VET18	Neostigmine 2.5 mg/mL	VET34	Hydromorphone 2 mg/mL
VET19	Lidocaine 20 mg/mL	VET35	Buprenorphine 0.3 mg/mL
VET20	Bupivacaine 5 mg/mL	VET36	Medathone ___ mg/mL
VET21	Furosemide 50 mg/mL	VET37	Norepinephrine ___ mcg/mL
VET22	Propofol 10 mg/mL	VET38	Dopamine 1600 mcg/mL
VET23	Alfaxalone 10 mg/mL	VET39	Dopamine 800 mcg/mL
VET24	Ketamine 100 mg/mL	VET40	Butorphanol 10 mg/mL
VET25	Atipamazole 5 mg/mL	VET41	Robenacoxib 20 mg/mL

## GENERIC DRUG LABELS

ITEM NO.	DESCRIPTION
VET42	Carprofen 50 mg/mL
VET43	Metoclopramide 5 mg/mL
VET44	Maropitant 10 mg/mL
VET45	Cefazolin 100 mg/mL
VET46	Meloxicam 20 mg/mL
VET47	Buprenorphine 0.5 mg/mL
VET48	Bupivacaine 5 mg/mL
VET49	Ondansetron 2 mg/mL
VET50	Pantoprazole 4 mg/mL
VET51	Enrofloxacin 50 mg/mL
VET52	Prepaid

**Don't see what you're looking for?**

Check out our extensive online catalogue with thousands of products to choose from!

**Pharmasytems.com**