



Dispensette & Supporting Supplies



Digital Dispensette

1 - 10 ml

- Easy, fast, and mechanically accurate liquid bottle-top dispenser
- Automated dispensing method
- Accurate to 0.5% of overall volumes
- Made with non-reactive Teflon®, resistant to wear and damage
- Ideal for dispensing liquid medication, especially methadone
- Recirculation tube prevents waste
- Safety valve prevents clogging and leaking
- Cap keeps out air, and prevents bacteria
- ISO/GLP compliant with certificate of pre-calibration; annual calibration recommended (contact us for details)

Item #: DK107-30

Lead time: 1-2 days

Minimum order quantity: 1 unit

Order Online



PharmaSystems Inc. | Your Single Source for Pharmacy Supplies



Important Maintenance Notes:

Includes 1 year warranty, 100% guarantee against product defects (does not include user-handling errors).

Maintenance for this equipment is very important; especially when dispensing methadone. When each stock bottle is complete, we recommend that you circulate the dispensette with hot water. Methadone crystallizes when exposed to air; without proper cleaning, it can clog the machine over time resulting in incorrect calibration. Should this occur, do not force your Dispensette. If your Dispensette becomes clogged over time, please contact PharmaSystems to schedule a cleaning. Refer to the instruction manual for full cleaning & care Instructions.

Dispensette includes certificate of calibration valid for 1 year after purchase. Annual calibration is recommended to ensure product accuracy. This service is offered by PharmaSystems for a fee. Please do not attempt to calibrate the Dispensette unless otherwise directed by a PharmaSystems employee as this will void your certificate of calibration. [Learn more about calibration >](#)



Digital Dispensette Extended Filling Tube

Extended filling tube required for use with 1L, 2L and 4L bottles

Item #: DK-201

Lead time: 1-2 days

Minimum order quantity: 1 unit

[Order Online](#)



1L Amber Glass Bottles

Recommendation to use at least 3 x 1L Amber Glass Bottles:

1. 1 for stock bottle
2. 1 to prepare next stock bottle
3. 1 for cleaning

Quantities: available as single units or in packs of 12

1 Bottle (item # 2001000-01)

12 per case (item # 2001000)

Lead time: 1-2 days

Order Online



White Lined Caps for use on 1L Amber Glass Bottles

- White caps, 33/400 neck finish.
- For use on 2001000 series 1L glass amber bottles

Item #: 20026

Pack of 25 caps

Lead time: 1-2 days

Order Online



Class A Glass Cylinders

- Meets ASTM E-1272 Class A, Style I Specifications & USP Standards for Volumetric Glassware
- Hexagonal base with 3 dimples for added stability - prevents cylinder from moving, even on wet surfaces
- Includes protective polyethylene hexagonal bumper guard
- Made of borosilicate glass
- Durable, white enamel graduations (applied by laser, heat annealed and inspected with code number to ensure accuracy)
- Volume accuracy calibrated to contain (TC)
- Recommended to be cleaned with distilled water to achieve good laboratory practice (GLP)

Item #: 1634AA-100

Lead time: 1-2 days

Minimum order quantity: 1 unit

Order Online



Cups

3.5oz Clear Plastic Cups

Item #: 7360

Quantity: packs of 100

Lead time: 1-2 days

[Order Online](#)

5oz Clear Plastic Cups

Item #: 7361

Quantity: packs of 50

Lead time: 1-2 days

[Order Online](#)



Warning Labels

Methadone Caution Label

Item #: P-2

Size: 1-9/16" x 3/8"

(40 x 9.5 mm)

Lead time: 1-2 days

2 rolls of 500 per box

Methadone may cause serious harm to someone other than the intended patient. Not to be used by anyone other than the patient for whom it was intended. MAY BE FATAL TO CHILD OR ADULT

[Order Online](#)

Methadone Caution Label (French)

Item #: FP-2

Size: 64 mm x 19 mm

Lead time: 1-2 days

Roll of 1,000

La méthadone peut sérieusement nuire à toute personne autre que le patient visé. Elle ne peut pas être utilisée par toute autre personne que le patient visé. PEUT ÊTRE FATALE À L'ENFANT OU À L'ADULTE

[Order Online](#)

Methadone can be fatal when taken by individuals for whom it is not prescribed...

PharmaSystems Item #:

CL-80

Size: 1-3/4 x 1" (44 x 25 mm)

Lead time: 1-2 days

2 rolls of 400 per box

Methadone can be fatal when taken by individuals for whom it is not prescribed. Shake well before using. May cause drowsiness. Keep medication stored in a locked box in the refrigerator. It is advised not to take alcoholic beverages when taking this medication.

[Order Online](#)

For Methadone Use Only

Item #: 13-2

Size: 44mm x 25mm

Lead time: 1-2 days

2 rolls of 400 per box

**FOR METHADONE
USE ONLY**

[Order Online](#)

Pricing is in CAD, FOB Markham, Ontario. Lead time is subject to change.



Warning Labels

Methadone 10 mg / ml

PharmaSystems Item #: CL-1150

Size: 1-9/16" x 3/8" (40 x 9.5 mm)

Lead time: 1-2 days

Pricing: \$42.00

(1,000 labels per box)

**Methadone
10 mg/ml**

[Order Online](#)

Methadone Label

PharmaSystems Item #: CL-81

Size: 1-9/16" x 3/8" (40 x 9.5 mm)

Lead time: 1-2 days

Pricing: \$42.00

(2 rolls of 500 per box)

METHADONE

[Order Online](#)

Pricing is in CAD, FOB Markham, Ontario. Lead time is subject to change.





View all dispensing supplies online [here >](#)

Connect with our bilingual customer service team for any questions or help placing your order:

www.PharmaSystems.com | info@pharmasystems.com

Phone: 905.475.2500 | Toll-Free: 1.888.475.2500

Fax: 905.475.7155 | Toll-Free: 1.888.475.7155

151 Telson Road, Markham Ontario Canada L3R 1E7

