



Globally Enhancing Patient Care

# Enteral Feeding Devices





---

# Contents

---

<b>Who we are</b>	<b>4</b>
<b>Enteral Syringe single use</b>	<b>6</b>
<b>Amber Barrel Syringes</b>	<b>8</b>
<b>Drainage and Venting</b>	<b>10</b>
<b>Stepped Bottle Adaptor</b>	<b>12</b>
<b>ENFit® Bottle Adaptor</b>	<b>13</b>
<b>Infant Feeding &amp; Breastfeeding Accessories</b>	<b>14</b>
<b>Tamper Evident Labels</b>	<b>14</b>
<b>ENFit® Bottle Transfer Cap</b>	<b>15</b>
<b>Enteral Caps</b>	<b>16</b>
<b>Push Fit Caps</b>	<b>17</b>
<b>Fixation Devices</b>	<b>18</b>
<b>Transition Adaptor</b>	<b>19</b>
<b>Enteral Connector Brush</b>	<b>20</b>
<b>pH</b>	<b>21</b>



**Distributed by PharmaSystems Inc.**

151 Telson Road, Markham, Ontario, Canada L3R 1E7

Please refer to the PharmaSystems product codes for information and pricing.  
Please note all product codes in grey are available by special order only.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)  
Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

---

# Who we are

---

**GBUK specialises in the supply of enteral feeding devices. Our comprehensive range of products enables healthcare professionals to safely deliver nutrition to their patients being cared for both in hospital and in their own homes.**

GBUK is proud to be playing a pivotal part in the global ISO ENFit® 80369-3 roll-out. We take pride in our commitment to quality, clinical innovation and providing exceptional value.

Committed to environmental sustainability with self generated on-site electricity from our wind turbine and solar panels.

## Our Mission:

GBUK, all subsidiaries and divisions are committed to improving patient outcomes and quality of life by delivering innovative, reliable, and easily accessible medical devices that empower healthcare professionals and serve global communities. As a trusted industry leader, we are committed to upholding the highest standards of regulatory compliance, ensuring safety and reliability in every product we offer, while creating lasting value for our customers and partners worldwide.

---

# Enteral Syringe

---

## ENFit® Syringe single use

ENFit® Single Use Syringes are used for the delivery of liquid medication, enteral feed and hydration.

Low dose connectors are available.



### Features and benefits

- Single use
- Clear barrel markings
- Full thread
- 1ml, 3ml and 5ml with low dose tip
- 60ml with eccentric tip
- 100ml syringe available
- DEHP and latex free

---

Distributed by PharmaSystems Inc.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)


Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

---

# Enteral Syringe

---

## ENFit® Syringe single use




<b>1ml</b>	500 pcs / box
Product Code:	<b>19798-31</b>



<b>3ml</b>	500 pcs / box
Product Code:	<b>19798-31</b>



<b>5ml</b>	500 pcs / box
Product Code:	<b>19800-31</b>



<b>10ml</b>	500 pcs / box
Product Code:	<b>19799-31</b>

---

Distributed by PharmaSystems Inc.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

---

# Enteral Syringe

---

## Amber Barrel

ISOSAF Amber barrel single-use syringes are compatible with most fluids including light sensitive or reactive, transmitting no more than 10% of incident light with wavelengths of 290–450nm.



## Features and benefits

- Single use
- Bulk non-sterile
- Includes push fit cap
- Clear barrel markings
- Full thread
- 1ml and 3ml with low dose tip
- DEHP and latex free



Includes  
Push Fit Cap




<b>1ml</b>	500 pcs / box
Product Code:	<b>20306-50</b>




<b>3ml</b>	500 pcs / box
Product Code:	<b>20307-50</b>




<b>5ml</b>	500 pcs / box
Product Code:	<b>20308-50</b>




<b>10ml</b>	500 pcs / box
Product Code:	<b>20309-50</b>



<b>20ml</b>	200 pcs / box
Product Code:	<b>20310-50</b>



<b>35ml</b>	100 pcs / box
Product Code:	<b>20311</b>



<b>60ml</b>	100 pcs / box
Product Code:	<b>20318</b>

**Distributed by PharmaSystems Inc.**


View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

# Drainage and Venting

## Enteral Drainage Bags

Enteral Drainage Bags fit any ENFit® compatible tubes and are used for gastric drainage.



## Features and benefits

- DEHP and latex free
- Male ENFit® and funnel connector available
- ENFit® compatible

EVA Drainage Bag			20 pcs / box Product Code:
50ml	Standard tubing	Distal ENFit connector	<b>DBISO-50</b>
100ml	Standard tubing	Distal ENFit connector	<b>DBISO-100</b>
250ml	Standard tubing	Distal ENFit connector	<b>DBISO-250</b>
500ml	Standard tubing	Distal ENFit connector	<b>DBISON-500</b>
	50cm tubing with clamp	Funnel connector	<b>DBISO-500</b>
1000ml	50cm tubing with clamp	Funnel connector	<b>DBISO-1000</b>

PVC Drainage Bag			50 pcs / box Product Code:
250ml	Standard tubing	Distal ENFit connector	<b>DBP-250</b>
500ml	Standard tubing	Distal ENFit connector	<b>DBP-500</b>
1000ml	Standard tubing	Funnel connector	<b>DBP-1000</b>
1500ml	Standard tubing	Funnel connector	<b>DBP-1500</b>
2000ml	Standard tubing	Funnel connector	<b>DBP-2000</b>


PVC Drainage Bag with two hooks			20 pcs / box Product Code:
250ml	50cm tubing with clamp	Distal ENFit connector	<b>DBP-250HT</b>
500ml	50cm tubing with clamp	Distal ENFit connector	<b>DBP-500HT</b>
1000ml	50cm tubing with clamp	Distal ENFit connector	<b>DBP-1000HT</b>

Distributed by PharmaSystems Inc.


View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

# Drainage and Venting



PVC Drainage Bags



PVC Drainage Bags with hooks



EVA Drainage Bags

Distributed by PharmaSystems Inc.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500


---

# Bottle Adaptor

---

## Stepped Bottle Adaptor

A dual-purpose, cost-effective medicine bottle adaptor that fits both ENFit® and oral syringes, drawing out the bottle contents.



### Features and benefits

- A cost-effective solution to last the lifetime of the bottle, once the bottle is empty discard the adaptor
- Polyurethane material moulds tightly to the bottleneck, for an improved fit
- Dual purpose – fits both ENFit® and oral syringes
- Tethered caps for secure storage
- Fits a wide range of bottle sizes which reduces the need to have different bottle adaptors

---

	50 pcs / box
Product Code:	<b>21052</b>

---


● **Tethered caps  
for secure storage**

● **Fits a wide range  
of bottle sizes  
which reduces the need  
to have different bottle  
adaptors**

# ENFit® Bottle Adaptor

## ENFit® Bottle Adaptor

ENFit® bottle adaptors are available in sizes to fit most medicine bottles, improving the dosage accuracy and reducing waste and contamination risk.



### Features and benefits

- Available in a range of sizes to fit the majority of medicine bottles
- Bottle cap can be replaced over adaptor in between uses
- Choice of packaging in bulk or single blisters

### Integrated

### tethered cap

the bottle adaptor can be left in place until the medicine is finished

Size:		100 pcs / box Product Code:
8-10mm	Single packed, non-sterile	<b>BAIS00-S</b>
12-15mm	Single packed, non-sterile	<b>BAIS01-S</b>
	Bulk non-sterile	<b>20312-10</b>
14-16.5ml	Single packed, non-sterile	<b>BAIS02-S</b>
	Bulk non-sterile	<b>20313-10</b>
17-20.5mm	Single packed, non-sterile	<b>BAIS03-S</b>
	Bulk non-sterile	<b>20314-10</b>
19-23mm	Single packed, non-sterile	<b>BAIS035-S</b>
	Bulk non-sterile	<b>20315-10</b>
24-26mm	Single packed, non-sterile	<b>BAIS04-S</b>
	Bulk non-sterile	<b>20316-10</b>


Distributed by PharmaSystems Inc.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

# Infant Feeding & Breastfeeding Accessories

## Tamper Evident Labels



### Features and benefits

- Anti-tamper scored edges that fray when peeled, meaning the label is unable to be removed from a baby bottle without clear signs of tampering
- The label also contains space to document patient and milk information
- Available in purple and yellow to distinguish between different milk types



35 pcs / box  
Product Code:

Purple

TEL-1

35 pcs / box  
Product Code:

Yellow

TELY-1

---

# Transfer Cap

---

## ENFit® Bottle Transfer Cap

An ENFit® compatible bottle transfer cap that enables the withdrawal of bottle content to a syringe, preventing leakage of contents.



### Features and benefits

- Tethered lid to protect contents, and for safe storage
- Air circulation feature, to avoid leakage of precious contents

### ● Closed system

A closed system to efficiently draw out the contents of the bottle using an ENFit® syringe

---

Product Code:	50 pcs / box
	<b>10180</b>

---

---

Distributed by PharmaSystems Inc.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

---

# Enteral Caps

---

## Male ENFit® Caps

ENFit® secure syringe caps with a male connector, allows for feed or medication to be made up in advance.



<b>Single Blister</b>	100 pcs / box
Product Code:	<b>ECISO</b>



<b>10 Caps / Pack</b>	100 pcs / box
Product Code:	<b>ECISO-10</b>

## Features and benefits

- Secure syringe caps to allow medication or feed to be prepared in advance
- Available in sterile packs of 10 or in single blisters

---

Distributed by PharmaSystems Inc.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500


---

# Push Fit Cap

---

## Enteral and Oral Push Fit Cap

ENFit® push fit cap that is compatible with ENFit® or oral syringes.



<b>Blister Sterile</b>	100 pcs / box
Product Code:	<b>PFC-100</b>
<b>Bulk Non-Sterile</b>	500 pcs / box
Product Code:	<b>20326</b>

### Features and benefits

- Weighted auto-stand cap will not fall connector side down if dropped
- Allows medication to be prepared in advance
- Allows for one handed use
- Self righting caps, allows for one handed capping
- Allows syringes to be capped off until required

<b>Bulk Non-Sterile</b>	200 pcs / box
Product Code:	<b>20325</b>

---

Distributed by PharmaSystems Inc.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

---

# Fixation Devices

---

## Clear Grip and Mini Clear Grip

Unique skin fixation devices to suit all tubes, ensuring secure and reliable placement.



### Features and benefits

- Ergonomic design ensures a comfortable fit, close to nostril to minimise tube pull outs
- Secures tubing between 2 layers of material, away from the patient's skin
- Clear material allows observation of skin

<b>Clear Grip Fixation Plaster</b>	30 pcs / box
--	--------------

Product Code:	<b>CG-30</b>
---------------	--------------

<b>Mini Clear Grip Fixation Plaster</b>	50 pcs / box
---	--------------

Product Code:	<b>MCG-50</b>
---------------	---------------



### Mini Clear Grip

Clear grip Mini has a silicone adhesive to protect the most delicate skin

---

Distributed by PharmaSystems Inc.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

---

# Transition Adaptor

---

## Transition Adaptor

Transition adaptors are ENFit® compatible, enabling connections during the transition to ENFit® ISO 80369-3.



### Features and benefits

- ENFit® compatible
- Latex free
- Tethered cap

● Cap enables closure of the tube in between access

100 pcs / box

Product Code:

**UAD2ISO**

---

Distributed by PharmaSystems Inc.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

---

# Connector Brush

---

## Enteral Connector Brush

The enteral connector brush will keep ENFit® connectors and caps free from debris and potential connection issues.



### Features and benefits

- ENFit® connectors have larger and deeper grooves, meaning increased area for debris to collect and potentially cause connection issues
- Using a brush after every use helps to prevent the accumulation of debris, and the feeding tube cap from adhering to the connector
- Keeping the connector and cap free of debris can help to extend the lifespan of the tube

In line with GEDSA ENFit cleaning Procedures 2018 recommendations Ref <http://stayconnected.org/wpcontent/uploads/2018/05/ENFit-Cleaning-Procedures-All-Tubes.pdf>

50 units / box

Product Code:

**ECB-NS**

---

# pH Testing

---

## Indicator strips 2-9

Enteral pH strips compatible with pHX-act pH reader



### Features and benefits

- Tethered cap
- Bright colour scale
- CE marked for the determination of human gastric aspirate
- Increments of 0.5
- Range pH 2-9
- 4 year shelf life
- Can be stored from 4°C to 30°C

### Patient Safety Alert

- In line with Patient Safety Alert NPSA/2011 PSA002 - Our pH indicator paper is CE marked and manufactured to test human gastric aspirate
- Supporting “Paediatric Nasogastric Tube Placement and Verification: Best Practice Recommendations from the NOVEL Project”

Irving, S. Y., Rempel, G., Lyman, B., Sevilla, W. M., Northington, L. Guenter, P. and , (2018), Paediatric Nasogastric Tube Placement and Verification: Best Practice Recommendations From the NOVEL Project. Nutrition in Clinical Practice, 33: 921-927. doi:10.1002/ncp.10189

	100 pcs / box
Product Code:	<b>ENT-PHX</b>

---

Distributed by PharmaSystems Inc.

View Pricing & Order Online: [www.PharmaSystems.com](http://www.PharmaSystems.com)

Email: [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com) | Call: 905.475.2500 | Toll Free: 1.888.475.2500

---

# pH Testing

---

## Indicator strips 2-9

A simple and effortless method for ascertaining the pH value of solutions with a colour scale you are familiar with.



### Features and benefits

- Tethered cap
- High quality pH paper
- CE marked for the determination of human gastric aspirate
- Increments of 0.5
- Accurate to pH +/-0-2

100 pcs / box

Product Code: **PHX2-9**

### Patient safety alert

- In line with Patient Safety Alert NPSA/2011 PSA002 - Our pH indicator paper is CE marked and manufactured to test human gastric aspirate
- Supporting “Paediatric Nasogastric Tube Placement and Verification: Best Practice Recommendations from the NOVEL Project”

Irving, S. Y., Rempel, G., Lyman, B., Sevilla, W. M., Northington, L. Guenter, P. and , (2018), Paediatric Nasogastric Tube Placement and Verification: Best Practice Recommendations From the NOVEL Project. Nutrition in Clinical Practice, 33: 921-927. doi:10.1002/ncp.10189





**GBUK Group Ltd.**

Woodland House, Blackwood Hall  
Business Park,  
North Duffield, Selby, North Yorkshire,  
YO8 5DD, UK

T +44 (0) 1757 288 587  
E [info@gbukgroup.com](mailto:info@gbukgroup.com)  
W [www.gbukgroup.com](http://www.gbukgroup.com)

**Distributed by PharmaSystems Inc.**

151 Telson Road  
Markham, Ontario  
Canada  
L3R 1E7

T 1.888.475.2500  
E [Info@PharmaSystems.com](mailto:Info@PharmaSystems.com)  
W [www.PharmaSystems.com](http://www.PharmaSystems.com)

