

Arisure®

Closed Male Luer (YM007) With Cap (YM104)

The Arisure Closed Male Luer (CML) is a safe and easy-to-use valved connector used for needle-free transfer of potentially hazardous fluids and medications between IV sets, syringes, and vials. When used with an Arisure Neutral Valve or other compatible valves, the Arisure CML forms an intuitive, secure Luer lock connection that enables a safe, simple transfer of medications¹. The Arisure CML enables a dry disconnection that reduces exposure to harmful drug residue. The Arisure CML remains mechanically and microbiologically closed, preventing leaks and reducing the risk of contamination.

¹ Data on file at Yukon Medical, LLC. ONE-LINK and CLEARLINK are registered trademarks of Baxter International, Inc. its subsidiaries, or affiliates. BD, the BD Logo, and SmartSite are trademarks of Becton, Dickinson and Company.

PRODUCT OVERVIEW

PRODUCT DESCRIPTION AND FUNCTION

The Arisure Closed Male Luer (CML) is single-use, sterile, non-pyrogenic, swab-able, normally closed, Luer-activated, bi-directional valved connector. It is used to transfer IV fluids when reconstituting, dispensing, administering, or disposing of potentially hazardous fluids or other medications.

The Arisure CML is versatile and can be used as a stand-alone male Luer connector on syringes or IV sets or with the Arisure Neutral Valve (and other compatible valves) during the handling of hazardous drugs to create a system that minimizes fluid leakage and the spread of vapor into the working environment.

INTENDED USE

Arisure is a Closed System Drug Transfer Device (CSTD) that mechanically prohibits the release of drugs in vapor, aerosol or liquid form during preparation and administration, and prevents the introduction of microbial and airborne contaminants into the drug or fluid path, allowing the system to minimize exposure of individuals, healthcare personnel, and the environment to hazardous drugs.

REGULATORY

US FDA Product Classification	ONB and FPA
CE marking	CE 2797
Quality management system certified	ISO 13485

PRODUCT CONFIGURATIONS

DEVICE CODE	DESCRIPTION	AVAILABLE FOR SALE	ORDERING INFO	PACKAGING INFO
 YM007	Arisure Closed Male Luer	US, CA, EU	Case Qty: 100	Case dimensions: 12in x 7.75in x 4in Gross weight per case: 1.6 lbs
 YM104 PharmaSystems Part Code: 88104	Arisure Closed Male Luer with Cap	US, CA	Case Qty: 100	Case dimensions: 12in x 7.75in x 4in Gross weight per case: 1.8 lbs