

NEW PHARMACY OPENING

2026 SUPPLY RECOMMENDATIONS



Quality Supplies and Exceptional Service

Proudly Canadian owned and operated since 1977, PharmaSystems is your Single Source for Pharmacy Supplies.

We have a Medical Device Establishment License for Class 1 and Class 2 devices, registered with Health Canada, and we are ISO Certified –13485:2016 MDSAP. We proudly support a wide range of healthcare & institutional partners across Canada.

We have created this curated catalogue to highlight our most-requested supplies to support new pharmacies opening.

Visit our full catalogue online:
PharmaSystems.com



Contents

4 Scales & Calibration

6 Weighing Accessories

7 Cylinders

10 Beakers, Stirring Rods

12 Flasks, Filter Paper

14 Spatulas

15 Mortar & Pestles

16 Ointment Accessories

18 Cups, Bottles & Vials

21 Prescription Bags, Tablet Counters

22 Prescription Management

25 Labels

26 Anti-Fatigue Mats

27 Thermometers

28 Order Form

How to order:

Online: www.pharmasystems.com
(register your account for access)

Call: 905.475.2500 (Toll free: 1.888.475.2500)

Fax: 905.475.7155 (Toll free: 1.888.475.7155)

Email: orders@pharmasystems.com

Our bilingual customer service team is readily available to help!

Office Hours: Monday-Friday
8:30am-1:00pm, 1:30-4:30pm, EST.

PharmaSystems Inc. 151 Telson Road
Markham Ontario Canada L3R 1E7

Shipping Terms: Orders will be shipped to all destinations in Canada. We reserve the privilege of shipping via the carrier of your choice.


Prices: Prices are subject to change without notice. Login to our website portal to view all catalogue pricing.

Returns: All returns to PharmaSystems must be authorized, and subject to a 25% restocking fee.

Accessibility: In meeting the needs of our customers and employees we strive at all times to provide our products in a way that respects the dignity of all individuals including persons with disabilities. In order to ensure equal treatment all PharmaSystems' customer service representatives and administrators are trained on policies, practices and procedures with respect to our commitment to people with disabilities.

Online Account Benefits:

- Competitive pricing
- 1000's of in-stock supplies ready to ship
- Quick inventory replenishment & same-day order processing
- Consolidate orders through a single source for supplies
- Shipping from our local warehouse, delivery anywhere in Canada



Breeze Break draft shield prevents air drifts caused by surrounding environmental movements - used to retain a constant mass reading!



EJ Series Electronic Scales

- Select from 3 decimal place display (VIC-123, VIC-303) or 2 decimal display (EJ200, EJ410, EJ610)
- Measures by weight, percentage, counting, and density functions
- Sample key available for easy parts counting, percentage and totals
- Outputs weight data, GLP report and calibration report to printer or computer (additional accessories are required, contact us for details)
- Full digital calibration – calibration weight not included (See item #CAL100/CAL200)
- Easy levelling – adjust 4 feet

Specifications	
Weighing Pan	110 mm (round)
Sensitivity Drift	± 20 ppm/ °C (10-30°C)
Min. % Display	0.1%
Display	7 segment LCD with backlight
Display Update	10 times per second
Operating Temperature	-10°C - 40°C, less than 85% R.H.
Power Supply	AC adapter (included) or 4 x "AA" size dry-cell batteries (70 hours of battery life)
Weight	850 g
Warranty	5 years

Item #	Description	Weight Capacity	Minimum Display	Accuracy	Weighing Pan Size	UPC
VIC-123 A	Compact Balance with Lid	120 g	0.001 g	±0.001 g	110 mm	063636112301
EJ200 B	Compact Precision Balance	210 g	0.01 g	±0.01 g	110 mm	063636510206
EJ200C	Compact Precision Balance with Lid	210 g	0.01 g	±0.01 g	110 mm	063636510206
VIC-303	Compact High Precision Balance with Lid	310 g	0.001 g	±0.001 g	110 mm	063636703035
EJ410	Compact Balance	410 g	0.01 g	±0.01 g	110 mm	063636514105
EJ610	Compact Precision Balance	610 g	0.01 g	±0.01 g	110 mm	063636516109
EJ-11 C	Breeze Break for use with EJ Series Balances	-	-	-	-	063636001100

A&D Precision Top Loading Balance

- Fast stabilization speed (1 second)
- Large bright vacuum florescent display
- Shock absorber function
- Statistical calculation function multiple weighing units
- Digital calibration function
- Counting function with ACAI (Automatic Counting Accuracy Improvement)
- Operating temperature - 5°C to 40°C (41°F to 104°F), 85% RH or less (No condensation)

Item #	Weight Capacity	Minimum Display	Accuracy	Weighing Pan Size	UPC
FX-300i A	320g	0.001g	±0.001g	130 mm	063636003005
FX-3000i B	3200g	0.01g	±0.001g	150 mm	063636783006
FZ-300I	320g	0.001g	±0.001g	130 mm	063636103002



Calibration Weights

- For calibration of electronic balance scales
- Recommended for use with scales to meet the requirements of the College of Pharmacists
- **Contact us for details on additional weights and scales that meet ASTM standards**

Item #	Description	UPC
CAL100 C	100g Calibration Weight	063636821104
CAL100GS	100g with Certificate of Accuracy	063636100124
CAL200	200g Calibration Weight	063636820008
CAL200GS	200g with Certificate of Accuracy	063636100131



Apothecary Weight Set

- Uniform density
- Low magnetic properties
- High resistance to wear and corrosion
- Includes: 10, 20, 50, 100, 200, 500 mg aluminum weights; ½, 1, 2, 3, 4, 5, grain aluminum links; 1, 2, 5, 10, 50 g.; 1, 2 oz., ½, 1 and dram, ½, 1 and 2 scruple brass weights; and tweezers


Item #	UPC
4027 D	063636740276



Octagonal Stir Bars

- PTFE coated
- Features ALNICO V magnetic cores (Store in pairs to maintain magnetic strength, away from scales)

Item #	Size	Length	Diameter	UPC
86301	Small	1"	5/16"	063636863012
86315	Medium	1.5"	5/16"	063636863159
86302	Large	2"	5/16"	063636863029



Paper Glassine Polished

- Non-absorbent glassine weighing paper
- 500 per package

Item #	Size	UPC
3410	3" square	063636
3016	4" square	063636301606
3412	6" square	063636
68415	6" square	063636784157

Plastic Weighing Boats

- Inexpensive, disposable containers
- Ideal for weighing quantities of ingredients
- Resilient, polystyrene construction ensures accurate control of pour-outs after weighing
- Non-stick surface is rugged enough for mixing
- Designed to easily stack and nest

Item #	Size	Dimensions	Qty	UPC
3420	Micro	4.4 x 4.4 cm	100	063636734206
40100	Small	4.5 x 4.5 cm	100	063636801007
3090	Medium	8.3 x 8.3 cm	100	063636730901
3092	Large	14 x 14 cm	100	063636730925



Anti-Static Weighing Canoes

- Provides anti-static surface for weighing ingredients
- Unique canoe shape enables easy pouring
- Made of virgin polystyrene with 10 to 9 OHM/square resistivity, the ideal dissipative range (the property of anti-static)
- Non-stick surface is rugged enough for mixing
- 100 per package

Item #	Size	Dimensions (inches)	UPC
3421	Small	1.8 x 2.2 x 5	063636734213
3422	Medium	3.25 x 5.28 x 1	063636734220
3423	Large	4.25 x 7 x 1.17	063636734237



Class A Graduated Cylinders


- High-quality autoclavable borosilicate glass cylinders
- Hexagonal base with 3 dimples for added stability
- White enamel graduations applied by laser, heat annealed, and inspected with code number to ensure accuracy
- Volume accuracy calibrated to contain (TC)
- Includes polyethylene hexagonal bumper guard (with the exception of the 5mL & 2000 mL size)
- Meets ASTM E-1272 Class A, Style I Specifications & USP Standards for Volumetric Glassware
- *Recommended for methadone dispensing

Item #	Volume	Increments	Accuracy	Scales	UPC
1634AA-5	5 mL	0.1 mL	±0.05 mL	Single	063636163426
1634AA-10	10 mL	0.1 mL	±0.05 mL	Single	063636163433
1634AA-25	25 mL	0.2 mL	±0.15 mL	Single	063636163457
1634AA-50	50 mL	0.5 mL	±0.25 mL	Double	063636163464
1634AA-100	A 100 mL	1 mL	±0.4 mL	Double	063636163488
1634AA-250*	250 mL	2 mL	±1 mL	Double	063636163495
1634AA-500	500 mL	5 mL	±2 mL	Double	063636163525
1634AA-1000	1000 mL	10 mL	±3 mL	Double	063636163532
1634AA-2000	2000 mL	20 mL	±6 mL	Double	063636163549

Class B Graduated Cylinders

- High-quality autoclavable borosilicate glass cylinders
- Detachable hexagonal sturdy base absorbs impact
- Includes protective hexagonal bumper guards
- White enamel graduations, applied by machine, heat annealed
- Volume accuracy calibrated to contain (TC)
- Meets ASTM E-1272 Class B, Style I Specifications

Item #	Volume	Increments	Accuracy	Scales	UPC
CYGPB-10	10 mL	0.1 mL	±0.1 mL	Single	063636882211
CYGPB-25	25 mL	0.2 mL	±0.3 mL	Single	063636881252
CYGPB-50	50 mL	0.4 mL	±0.4 mL	Double	063636818500
CYGPB-100	B 100 mL	1 mL	±0.6 mL	Double	063636817107





Bell & Conical Glass Cylinders

- Bell and conical shaped glass graduated cylinders, with easy-pour spouts
- Bell shape has a controlled diameter – recommend for accuracy
- Bell-shaped cylinders have single scale white enamel graduations, manually applied and annealed for accuracy
- Conical-shaped cylinders have single scale molded graduations
- Volume accuracy calibrated to contain (TC)
- Flint conical glass cylinders cannot be autoclaved

Item #	Volume	Shape	Material	Increments	UPC (063636)
3360	5 mL	Bell	Borosilicate	1 mL	063636833602
3361	10 mL	Bell	Borosilicate	0.5 mL	063636833619
3362	25 mL	Bell	Borosilicate	1 mL	063636833626
3019	50 mL	Conical	Flint	5 mL	063636830199
3021	100 mL	Conical	Flint	5 – 25 mL	063636830212
3022	250 mL	Conical	Flint	10 – 50 mL	063636830229
3023	A 500 mL	Conical	Flint	25 – 50 mL	063636830236
3363	B 1000 mL	Bell	Borosilicate	20 mL	063636733636

Glass Cylinders

- High-quality autoclavable borosilicate glass
- White enamel graduations, applied by machine, heat annealed
- Volume accuracy calibrated to contain (TC)
- Round base
- Meets Class B specifications

Item #	Volume	Increments	Scales	UPC
3002	10 mL	0.2 mL	Single	063636830021
3009	25 mL	0.5 mL	Single	063636830090
3006	50 mL	1 mL	Single	063636830069
3008	100 mL	1 mL	Single	063636830083
3014	C 250 mL	2 mL	Single	063636830144
3015	500 mL	5 mL	Single	063636830151
3371	1000 mL	10 mL	Single	063636833718

Clear Plastic Graduated Cylinders

- Polymethylpentene (TPX) plastic cylinders are highly resistant to chemical attack (at room temperature)
- Sturdy hexagonal base
- Clear plastic with molded graduations
- Can be sterilized by autoclaving (120°C), gas, and radiation – not autoclavable by dry heat
- Can hold a maximum temperature of 175°C
- Volumes calibrated to deliver (TD)

Item #	Volume	Increments	Accuracy	Scales	UPC
3040	10 mL	0.2 mL	±0.2 mL	Single	063636830403
3435	25 mL	0.5 mL	±0.5 mL	Single	063636834357
3045	50 mL	1 mL	±1 mL	Single	063636830458
3046	100 mL	1 mL	±1 mL	Single	063636830465
3047	A 250 mL	2 mL	±2 mL	Single	063636830472
3048	500 mL	5 mL	±5 mL	Single	063636830489
3436	1000 mL	10 mL	±10 mL	Single	063636834364
3437	2000 mL	20 mL	±20 mL	Single	063636834371

Plastic Graduated Cylinders

- Polypropylene cylinders are very resistant to effects of Alkalies
- Sturdy hexagonal base
- Can be sterilized by autoclaving (120°C), gas, and radiation – not autoclavable by dry heat
- Can hold a maximum temperature of 135°C
- Volumes calibrated to deliver (TD)

Item #	Volume	Increments	Accuracy	Scales	UPC
3431	10 mL	0.2 mL	±0.2 mL	Single	063636834319
3034	50 mL	1 mL	±1 mL	Single	063636830342
3036	100 mL	1 mL	±1 mL	Single	063636830366
3037	B 250 mL	2 mL	±2 mL	Single	063636830373
3038	500 mL	5 mL	±5 mL	Single	063636830380
3439	2000 mL	20 mL	±20 mL	Single	063636343903

Cylinder Brushes

- Use to clean cylinders, graduates and bottles
- Sturdy wire handled brush with nylon bristles and tufted end for thorough cleaning
- Size includes handle

Item #	Description	Size	UPC
3049	C Small - For 50 mL Cylinders	11" (229 mm)	063636730499
3059	Medium - For 100-200 mL Cylinders	13" (330 mm)	063636730598
3066	Large - For 500 mL Cylinders	17" (438 mm)	063636730666



Pitchers

- Transparent polypropylene plastic
- Handles for convenient lifting
- Autoclavable and chemically resistant
- Raised graduations

Item #	Volume	Increments	UPC
PIP-500	500 mL	5 mL	063636747053
PIP-1L	1000 mL	10 mL	063636747107
PIP-2L	2000 mL	50 mL	063636747206
40055	3000 mL	100 mL	063636700553
PIP-5L	5000 mL	250 mL	063636747503



Clear Plastic Beakers

- Polymethylpentene (TPX) is highly resistant to chemical attack (at room temperature)
- Can be sterilized by autoclaving (120°C), gas, and radiation – not autoclavable by dry heat
- Molded graduations

Item #	Volume	Increments	UPC
BKTPX-50	50 mL	2 mL	063636790509
BKTPX-100	100 mL	5 mL	063636891008
BKTPX-250	250 mL	10 mL	063636792503
BKTPX-600	600 mL	10 mL	063636896003
BKTPX-1L	1000 mL	20 mL	063636740016
BKTPX-2L	2000 mL	50 mL	063636740023



Glass Beakers

- Low-form Griffin borosilicate glass beaker with white enamel single scale graduations, applied by laser and heat annealed
- Shape and heavy wall construction offers optimum strength
- Spout for easy pouring
- Ideal for mixing and heating materials

Item #	Volume	Increments	UPC
3372	30 mL	10 mL	063636833725
3024	50 mL	10 mL	063636830243
3025	100 mL	10 mL	063636830250
3373	150 mL	25 mL	063636833732
3026	250 mL	25 mL	063636830267
3374	400 mL	25 mL	063636833749
3027	600 mL	50 mL	063636830274
3375	1000 mL	50 mL	063636833756



Plastic Beakers

- Graduations are molded in milliliters
- Polypropylene beakers are very resistant to effects of Alkalies
- Can be sterilized by autoclaving (120°C), gas, and radiation – not autoclavable by dry heat

Item #	Volume	Increments	UPC
BKP-50	50 mL	2 mL	063636700508
BKP-100	100 mL	5 mL	063636821005
BKP-250	250 mL	10 mL	063636822507
BKP-600	600 mL	10 mL	063636826000
BKP-1L	1000 mL	20 mL	063636820015
BKP-2L	2000 mL	50 mL	063636820022
BKP-5L	5000 mL	500 mL	063636820053

Stainless Steel Beakers

- Durable, seamless construction, with easy-pour rim
- Will not break, leach product, or become marred inside
- Non-corrosive and rust resistant

Item	Description	UPC (063636)
18602	600 mL	186029
18604	2000 mL, with handle	186050
18603	4000 mL, with handle	186036



Glass Stirring Rods

- Soft glass with rounded ends minimizes scratching

Item	Size	Qty	UPC
3017	8"	3	063636830175
3018	10"	3	063636830182

Plastic Stirring Rods

- White polypropylene plastic
- Small and large spatulas form at each end

Item #	Size	UPC
39160	9.63" (24.5 cm) long	063636791605

Stir Bar Retriever

- Ideal for use with stirring rods (see above)
- Encased in chemically resistant polypropylene
- Hand ring allows for hanging when not in use

Item #	Description	UPC
86398	18 inches (45.72 cm) long	063636863982



Erlenmeyer Flasks

- White permanent graduations and large marking squares
- Made from high-quality borosilicate tubing with a low coefficient of expansion
- Volumes calibrated to deliver (TD)
- Meets ASTM E-1404 Type I, Class 1 specifications



Item #	Volume	Increments	UPC
24-125	125 mL	25 mL	063636741259
24-250	250 mL	50 mL	063636742508
24-500	500 mL	100 mL	063636745004
24-1L	1000 mL	100 mL	063636710002
24-2L	2000 mL	250 mL	063636720001
24-5L	5000 mL	500 mL	063636750008



Volumetric Flasks

- Class A volumetric flask, calibrated to contain (TC) with fitted polyethylene stopper
- Borosilicate tubing provides uniform wall thickness resulting in optimum balance between thermal shock resistance and mechanical strength
- White enamel graduations

Item #	Volume	Accuracy	UPC (063636)
2415-10	10 mL	±0.02 mL	063636241513
2415-25	25 mL	±0.03 mL	063636241520
2415-50	50 mL	±0.05 mL	063636241551
2415-100	100 mL	±0.1 mL	063636751005
2415-200	200 mL	±0.15 mL	063636224158
2415-250	250 mL	±0.15 mL	063636250249
2415-500	500 mL	±0.25 mL	063636755003
2415-1000	1000 mL	±0.4 mL	063636750015
2415-2000	2000 mL	±0.6 mL	063636222413

Plastic Funnels

- Lightweight polypropylene
- Durable, long stem design

Item #	Diameter	Capacity	UPC
3030	55 mm	55 mL	063636830304
3031	100 mm	100 mL	063636830311

Heavy Duty Glass Funnel

- Ribbed on inner surface to increase flow

Item #	Diameter	Capacity	UPC
3061	7 cm	60 mL	063636730611
3063	9 cm	125 mL	063636730635
3065	11½ cm	250 mL	063636730659
35695-04	14½ cm	500 mL	063636795047
35695-05	14½ cm	100 mL	063636795054



Filter Paper

- Use for filtration of coarse particles and gelatinous precipitates
- Standard qualitative filter papers suitable for quadrant fold or gravity flow applications
- Contains high level of alpha-cellulose that provides high purity, weakens when wet

Item #	Diameter	Grade	Particle Retention	Flow Rate	Speed	Rapidity	Qty	UPC
2994	110 mm	04	10 µm	Medium	34	200 mls/min	100	063636729943
2995	75 mm	04	10 µm	Medium	34	200 mls/min	100	063636729950
2998	110 mm	01	2.5 µm	Medium	124	55 mls/min	100	063636029999
2999	185 mm	01	2.5 µm	Medium	124	55 mls/min	100	063636029982
2996	240 mm	01	2.5 µm	Medium	124	55 mls/min	100	063636729967
2997	320 mm	01	2.5 µm	Medium	124	55 mls/min	100	063636729974

pH Test Papers

Item # 90710

- 1-14 pH range test papers
- Supplied in a clear plastic vial container
- Includes colour guides in vial
- 100 test strips per vial
- Each test strip measures 3 x 12 mm

Item # 6080

- 6-8 pH range
- Double roll, hydron pH test papers
- Includes double roll dispenser
- Includes pH colour charts

Item #	Description	UPC
90710	A 1-14 pH range, 100 per vial	063636707101
6080	B 6-8 (neutral) pH range, double roll with dispenser	063636760809

Instructions for use of the test strip (item # 90710):

1. Use first morning urine
2. After a steady stream has started, wet the strip pad for 1-2 seconds
3. Gently shake off any excess fluid on the test strip and wait for 10-15 seconds
4. Compare/match the color of the test strip to the colour chart provided in the bottle. The number that corresponds to the color matched pad is the result.

The normal pH of urine is slightly acidic and gives a reading between 6.5-7.5. If the urine is on the alkaline side, the readings should correspond to colors indicating 7.5 or higher. Acidic urine gives color corresponding to 6.5 or lower. Please note that pH of urine is affected by the food and medication intake.



Rubber Spatulas

- Hard, non-porous rubber spatula for compounding
- Recommended for use with compounds that react with metal
- Can be cleaned with soap and water, alcohol or by autoclaving

Item #	Blade Size	UPC
3091	A 4" (101.6 mm)	063636830915
3093	6" (152.4 mm)	063636830939
3095	8" (203.2 mm)	063636830953

Spatula with Multi-Purpose Tool

- Includes swing-out multi-purpose cutting tool
- Tapered end ideal for use when counting tablets
- Blade and cutting tool are made of stainless steel, handle is made of rosewood
- Cutting tool measures 1 1/4"

Item #	UPC
1137	B 063636711375

Stainless Steel Spatulas


- Flexible, dull curved edged, with a comfortable wood handle
- 4", 6" and 8" spatulas are recommended to meet the requirements of the College of Pharmacists

Item #	Blade Size	UPC
3042	4" (101.6 mm)	063636304249
3043	5" (114.3 mm)	063636304348
3028	C 6" (152.4 mm)	063636830281
3029	8" (203.2 mm)	063636830298
3429	10" (254 mm)	063636834296
3430	12" (355.6 mm)	063636834302

Disposable Spatulas

- Unique hollow cores minimize thermal transfer
- Resistant to dilute acids and bases
- Autoclavable polypropylene can be used for sterile compounding

Item #	Description	Size	Qty	UPC
17946	D Disposable, Opaque	210mm	300	063636179465
17947	E Disposable, Blue	210mm	300	063636179472
17948	F Disposable, Green	310mm	150	063636179489



Flint Glass



Wedgewood

Flint Glass Mortar & Pestles

- Convenient - Use for preparing medications for patients or compounding recipes.
- Easy to Clean - Can be sterilized, cleaned with alcohol or washed with soap and water.
- **Note:** Not recommended for use with drastic temperature changes.

Wedgewood Mortar & Pestles

- Convenient - Use for preparing medications for patients or compounding recipes.
- Easy to Clean - Can be cleaned with alcohol or washed with soap and water. Pestles cannot be sterilized because of the wooden handle.
- **Note:** Wedgewood should not be used with flavors or scents because the porous material will retain these elements after they are cleaned.

Mortar & Pestles

Item #	Description	Size	UPC
3070	Flint Glass Mortar	2 oz. (60 mL)	063636830700
3080	Flint Glass Pestle	2 oz. (60 mL)	063636830809
3071	Flint Glass Mortar & Pestle Set	2 oz. (60 mL)	063636830717
3072	Flint Glass Mortar	4 oz. (120 mL)	063636730727
3082	Flint Glass Pestle	4 oz. (120 mL)	063636730826
3073	Flint Glass Mortar & Pestle Set	4 oz. (120 mL)	063636830731
3074	Flint Glass Mortar	8 oz. (240 mL)	063636730741
3084	Flint Glass Pestle	8 oz. (240 mL)	063636730840
3075	Flint Glass Mortar & Pestle Set	8 oz. (240 mL)	063636830755
3076	Flint Glass Mortar	16 oz. (480 mL)	063636730762
3086	Flint Glass Pestle	16 oz. (480 mL)	063636830861
3382	Flint Glass Mortar & Pestle Set	16 oz. (480 mL)	063636033828
3053	Wedgewood Mortar & Pestle set	4 oz. (120 mL)	063636730536
3055	Wedgewood Mortar & Pestle set	8 oz. (240 mL)	063636730550
3395	Ceramic Mortar & Pestle set	2 oz. (60 mL)	063636733957
3398	Ceramic Mortar & Pestle set	4 oz. (120 mL)	063636733988
3401	Ceramic Mortar & Pestle set	8 oz. (240 mL)	063636734015
3404	Ceramic Mortar & Pestle set	16 oz. (480 mL)	063636734046
3407	Ceramic Mortar & Pestle set	32 oz. (960 mL)	063636734077

Ceramic Mortar & Pestles

- Use for preparing medications for patients or mixing and pulverizing ingredients for compounding
- Ceramic models are ideal for compounding tablets and solids
- Can be sterilized, cleaned with alcohol, or washed with soap and water



Ceramic

Ointment Pad


- 12" square
- 100 sheets per pad

Item #	Description	Qty	UPC
3088	A 12 x 12" Square	100 Sheets	063636308803

Glass Ointment Slab

- Glass ointment slab suitable for mixing various preparations
- 12" square
- 1/4" thick

Item #	Description	UPC
3098	B 12" Square Slab, 1/4" Thick	063636830984



Ointment Tubes

- Aluminum Ointment Tubes (A-F):
- Ideal for use with creams, ointments and topical medications
- White metal tube with white plastic screw cap and plastic ophthalmic tip - tubes lined with Teflon® to help prevent sticking or reactions with the aluminum; autoclavable
- Secure by folding the ends with tongs or pliers

Item #	Description	Qty	UPC
10205-01	C 5 g Aluminum Ointment Tube, with Tamper Evident Ring	135	063636702052
10206-01	D 10 g Aluminum Ointment Tube	100	063636702069
10201-01	E 15 g Aluminum Ointment Tube	117	063636702014
10202	F 30 g Aluminum Ointment Tube	12	063636702021
10203A	G 50 g Aluminum Ointment Tube	12	063636702038
10204-01	H 100 g Aluminum Ointment Tube	54	063636720414
4040	I Ointment Tube Pliers	1	063636740405

Dispensing Containers and Accessories

- Built-In Safety - Bottom closure snaps shut for a tamper-resistant seal
- Maintain Freshness - Screw top lids keep ointments fresh until the last application
- Easy to Use - Simply twist the bottom of the container and the compound is expelled from the top of the container

Item #	Description	Qty	UPC
10332	A 30 g Dispensing Container	20	063636103323
10333	50 g Dispensing Container	20	063636103330
10334	100 g Dispensing Container	15	063636103347
10335	B 200 g Dispensing Container	10	063636103354
10330	C Large Mixing Disc	30	063636103309



Ointment Jars

- Smart Packaging - Containers come shrink wrapped on paperboard trays to keep them clean and provide convenient storage
- Multi-Purpose - Ideal for packaging and dispensing ointments, gels and creams

Item #	Size	Qty	UPC
10083	25 mL	120	063636700836
10084	40 mL	80	063636700843
10085	55 mL	80	063636700850
10086	A 75 mL	35	063636700867
10087	120 mL	35	063636700874
10088	170 mL	35	063636700881
10089	B 250 mL	24	063636700898
10520	300 mL	10	063636705206
10521	C 690 mL	10	063636705213

Medication Cups

Item #	Description	Type	Qty	UPC
7364	1 oz. / 29 mL Graduated	Plastic	1000	063636873646
7362	2 oz. / 59 mL Graduated	Plastic	1000	063636127626
5165	1 oz. / 29 mL Graduated	Plastic	400	063636751654
5165-01	1 oz. / 29 mL Graduated	Plastic	4800	-
7360	3.5 oz. / 103 mL	Plastic	100	063636126612
7361	5 oz. / 147 mL	Plastic	100	063636127619
7366	5 oz. / 147 mL	Plastic	1200	063636773663
7369	1 oz / 29 mL	Paper Souffle	1000	063636873691
7368	1/2 oz / 14 mL	Paper Souffle	1000	063636873684
7363	3 oz / 88 mL	Wax Paper	2400	063636073633
7363	3 oz / 88 mL	Wax Paper	2400	063636073633
5160	Lids for Narrow Cups	Plastic	400	063636751609
5160-01	Lids for Narrow Cups	Plastic	4000	063636751609



Amber Glass Medicine Bottles

- 28 mm white plastic screw caps included

Item #	Description	Bottle Neck Size	Size	Qty	UPC
20050-12	Amber Glass	18 mm	50 mL	264	-
200100	Amber Glass	18 mm	100 mL	140	063636200121
200200	Amber Glass	24 mm	236 mL	108	-
200500-12	Amber Glass	28 mm	500 mL	12	-
2001000	Amber Glass	33 mm	1 L	12	063636200114
2001000-01	Amber Glass	33 mm	1 L	1	063636710170
20020	White Cap For use with 20050 & 200100			100	063636200213
20026	White Cap For use with 2001000 & 2001000-01			25	063636200206



Amber Plastic Bottles, PET

- Bottle necks measure 24mm – compatible with E-28mm and F-28mm Adapta-Caps (see page 49 for details)

Item #	Description	Includes Caps	Qty	UPC
102893	100 mL Amber Plastic Bottles	✓	108	063636728939
10289	100 mL Amber Plastic Bottles	-	108	063636702892
102903	200 mL Amber Plastic Bottles	✓	63	063636729035
10290	200 mL Amber Plastic Bottles	-	63	063636702908
102913	300 mL Amber Plastic Bottles	✓	48	063636729134
10291	300 mL Amber Plastic Bottles	-	48	063636702915
102923	500 mL Amber Plastic Bottles	✓	35	063636729233
10292	500 mL Amber Plastic Bottles	-	35	063636702922
102933	1000 mL Amber Plastic Bottles	✓	30	063636729332
10293	1000 mL Amber Plastic Bottles	-	30	063636702939
10294	Tamper Resistant Caps (Not child-proof)	-	250	063636702946



Amber Glass Dropper Bottles with Glass Dropper Caps

Item #	Description	Size	Qty.	UPC (063636)
19450	Amber Glass with Dropper Cap	5 mL	340	063636194505
19451	Amber Glass with Dropper Cap	10 mL	255	063636194512
19452	Amber Glass with Dropper Cap	20 mL	216	063636194529
19453	Amber Glass with Dropper Cap	30 mL	168	063636194536
19454	Amber Glass with Dropper Cap	50 mL	99	063636194543
19455	Amber Glass with Dropper Cap	100 mL	68	063636194550

Amber Plastic Dropper Bottles

- Bottle caps have built-in droppers for dispensing or mixing
- Material: Polypropylene

Item #	Size (oz.)	Size (mL)	Qty	UPC
10078	½ oz	14 mL	12	063636800789
10079	1 oz	29 mL	12	063636800796
10080	2 oz	59 mL	12	063636800802



Amber Glass Dropper Bottle Accessories

Item #	Description	Qty.	UPC
10067	A White Tamper Resistant Caps	100	063636100674
10068	Black Tamper Resistant Caps	100	063636100681
19456	B Glass Dropper or 5 mL Bottles	100	063636194567
19457	Glass Dropper for 10-30 mL Bottles	100	063636194574
19458	Glass Dropper for 50-100 mL Bottles	100	063636194581
19459	C Glass Dropper for item # 19453	100	063636194598



Screw-on caps have a tamper-evident ring that breaks off to show evidence of access.

Items highlighted in orange are recommended for methadone dispensing



Medicine Vials

- Translucent polyethylene amber vials provide visibility
- Compatible with child-resistant or snap caps (see options below)
- When closed with the Child Resistant Caps, vials pass both child resistant and senior friendly protocol tests
- Use with either caps to meet USP31/NF26 requirements
- Vials and caps pass USP671 Moisture Permeation tests

Item #	Description	Vial Size/Compatible with Vial	Qty	UPC
572532	Medicine Vial with Child Resistant Cap Attached	33 mL (9 dram)	100	063636725327
572633	Medicine Vial with Child Resistant Cap Attached	48 mL (13 dram)	100	063636726331
572733	Medicine Vial with Child Resistant Cap Attached	60 mL (16 dram)	100	063636727338
572834	Medicine Vial with Child Resistant Cap Attached	73 mL (20 dram)	100	063636828349
572935	Medicine Vial with Child Resistant Cap Attached	110 mL (30 dram)	100	063636729356
573035	Medicine Vial with Child Resistant Cap Attached	148 mL (40 dram)	100	063636730352
573135	Medicine Vial with Child Resistant Cap Attached	220 mL (60 dram)	100	063636731359
572536	Medicine Vial with Snap Cap Attached	33 mL (9 dram)	100	063636725365

Sterile Dropper Bottles

- Sterility Assurance - Ideal for the preparation of sterile ophthalmic solutions. Bottle exterior is sterile until the package is opened; the inside is sterile until the cap is removed.
- All-in-One System - Includes snap-in dropper tips and screw on caps.
- Controlled Dispensing - Releases approximately 42uL per drop and the dispensing tolerance is +/- 5uL

Item #	Size	Amount dispensed per drop	Qty	UPC
7780	3 mL	0.023 mL	12	063636777807
7780-02	3 mL	0.023 mL	2	063636877804
7782	7 mL	0.035 mL	12	063636777821
7782-01	7 mL	0.035 mL	1	063636877828
7784	15 mL	0.035 mL	12	063636777845
7784-01	15 mL	0.035 mL	1	063636877842
7786	30 mL	0.035 mL	12	063636777869
7786-01	30 mL	0.035 mL	1	063636877866



Prescription Bags

- Convenient - Printed white paper bags in three sizes provide a convenient, cost-effective method of med packaging at the pharmacy and upon patient discharge. Largest size bag has a flat bottom for convenient filling.
- No Additional Fees - No order minimums and no set-up fees.
- Package - 100 per package

Item #	Dimensions (W x H x D)	UPC
19927	A 3 1/2 x 9 3/4 x 1 1/2"	063636199272
19928	B 5 x 10 x 2"	063636199289
19929	5 1/4 x 10 15/16 x 3 7/8"	-



Fretwell's Triangle

- Count regular shaped tablets with absolute accuracy and speed
- Tablets glide easily on the stainless steel triangle
- Easy to clean

Item #	UPC
FT1	C 063636886202

Tablet Counters

- Plastic pill counting tray with clear lip handle and two corner spouts
- Effective way to keep organized with pill counting
- Includes plastic spatula (0.6 x 5 1/2") to aid counting and transferring pills
- Recommended for right-handed users
- Measures: 8 1/4 x 6 x 1"

Item #	Colour	UPC
D-50	D Clear/Blue	063636050306

Ergonomic Tablet Counters

- Includes spatula with ergonomic grip for easy handling and docking slot
- Designed to provide stability while counting tablets
- Deep blue colour allows pills and powder traces to contrast with tray
- Measures: 7 3/4 x 6.15 x 0.4"

Item #	Description	UPC
7147	Right-Handed	063636771478
7148	E Left-Handed	063636771485

Ambidextrous Tablet Counter

- Polyethylene pill counting tray for left or right handed users
- Spatula included
- Chute has a hinged cover for easy filling
- Measures: 7 3/4 x 6 1/4 x 1 1/2"

Item #	Description	UPC
12261	F Left or Right Handed	063636122614

Hanging Prescription Bags

- Simplify prescription organization and retrieval
- Made of tough, clear 4 mil plastic for extreme durability
- Fits on any rod with a diameter of 1 cm or smaller

Item #	Size (W x H)	Qty	UPC
17540 A	10 1/2 x 7 1/2"	10	063636175405
17541	12 x 9"	10	063636175412
17542	14 x 12 1/2"	10	063636175429

Wall Rack for Prescription Bags

- Simple solution to prescription organization
- Saves time, improves accuracy and enhances customer satisfaction
- Aluminum frame is 3 ft. wide and holds five 1 cm thick rods
- Each rod holds approximately 8-10 bags per foot or as many as 150 bags per rack
- Easy assembly, all hardware included
- Measures: 36" W x 57 1/2" H x 7 1/2" D
- Use with hanging bags and ID tags (sold separately)

Item #	Description	Qty	UPC
17543	Wall Rack for Prescription Bags	1 frame, 5 rods	063636175436
17544 B	Green ID Tags	25	063636175443



Stock Prescription Receipt Pads

- Three popular sizes and designs
- Padded in 100's, sold in 1000's

Item #	Description	Qty	UPC
PRP-3 A	4 1/4 x 7"	10 pads of 100 sheets	063636300005
PRP-6 B	3 1/4 x 7"	10 pads of 100 sheets	063636660000
PP-8 C	4 1/4 x 6" (Hole Punched)	10 pads of 100 sheets	063636800000



Rx Tape

- 3 sizes are individually packaged for convenient storage
- Two core sizes to fit different tape dispensers

Item #	Width	Length	Core Size	UPC
T1-1	1"	33 yards	1"	063636024338
T1-3	1"	72 yards	3"	063636024666
T1-5	1 1/2"	72 yards	3"	063636036669
T1-7	2"	55 yards	3"	063636048501



Multi-Cut Utility Shears

- Stainless steel blade, ideal for cutting bandages, dressing, fabric, or other tough material
- Colours available: ● black, ● blue, ● orange, ● pink, ● purple, ● red and ● yellow -specify when ordering

Item #	UPC
95-105 C	063636941956



Handle Universal Scissors

- Stainless steel blade

Item #	Size	UPC
S3115 D	5 1/2"	063636831158



Forceps

- Made of stainless steel

Item #	Description	Size	UPC
F4111	Dressing Forceps	4 1/2"	063636841119
3039	Dressing Forceps	5 1/2"	063636730390
F0111 E	Super Fine Straight Splinter Forceps	4 1/2"	063636701116

Tape Dispensers

Single Roll Tape Dispenser

- Accepts tape up to 1" wide
- Core adapter accommodates 1" and 3" cores

Double Roll Tape Dispenser

- Accepts tape up to 2" wide
- Core adapter accommodates 3" cores

Triple Roll Tape Dispenser

- Accepts tape up to 3" wide
- Can run 3 rolls of 1" wide tape or two rolls of 1 1/2" tape individually, 3" core

Item #	Description	UPC
TD1 D	Single Roll Tape Dispenser	063636010003
GF2306 E	Double Roll Tape Dispenser	063636230661
GF2336 F	Triple Roll Tape Dispenser	063636233662



Prescription Boxes

- Storage - Use for storing archives or keep it on the counter for organizing incoming prescriptions
- Shipping - The paperboard box is shipped flat and is easy to assemble.

Item #	Size L x W x H	Qty	UPC
11194 A	Small 13 x 6 ¹ / ₄ x 4 ¹ / ₂ "	1	063636711948
10196 B	Large 16 ¹ / ₈ x 5 ³ / ₄ x 7 ¹ / ₂ "	1	063636701963



Prescription File Folders

- Keeps prescription slips organized for quick reference and storage
- Double fold file folders hold approximately 100 prescription sheets
- When folded, the spine is 1/2" deep and can be used to label or identify contents
- Measures: 5 1/2 x 13"
- Accommodates: 5 1/2 x 4" prescription slips



Item #	Qty	UPC
1976 C	100	063636819767
19464 D	100	063636194642

Tongue Depressors

- Sterile
- Individually wrapped tongue depressors
- Sterilized and wrapped in peel open envelope
- 100 per pack

Item #	Qty	UPC
T1100	100	063636110000



Label Dispensers


- Rugged construction and high quality design ensure years of trouble free use
- Superior design permits one hand label dispensing
- Dispensers can be quickly loaded and easily cleaned
- Holds any of the PharmaSystems® standard size warning labels, shown on previous pages in this section
- Additional, pre-loaded options below



Item #	Description	UPC
BK-10	Label Dispenser for 10 Rolls	063636000103
BK-20	Label Dispenser for 20 Rolls	063636000202
BK-30	Label Dispenser for 30 Rolls	063636000301

17/20, Warning Label System

- 17 rolls of 500 auxiliary labels preloaded into a 20 roll label dispenser (includes BK-20 model)
- Select from English or French sets



Item #	Description	UPC Code
BK-40 A	English Warning Label System	063636404048
FW-1720 B	French Warning Label System	063636172022

BK-40	FW-1720
P-1	FG-61
1X	FW-15
2X	F-669
3W	FX-1
4X	FX-10
5X	FX-11
6X	FX-12
P-12	FX-2
8X	FX-5
1-10G	FX-4
10X	FG-27
11X	FX-8
12X	FG-10
1-15G	F-690
15W	FG-77
16W	FG-13
1-11G	FG-15

20/30, Warning Label System

- 20 rolls of 500 auxiliary labels preloaded into a 30 roll label dispenser (includes BK-30 model)
- Select from English or French sets



Item #	Description	UPC Code
BK-50 A	English Warning Label System	063636494940
FW-1730 B	French Warning Label System	063636173036

BK-50	FW-1730
1X	FG-15
2X	FW-15
3W	FW-3
4X	FX-1
5X	FX-10
8X	FX-11
10X	FX-12
11X	FX-2
15W	FX-4
1-10G	FX-5
1-11G	FG-61
1-49G	FX-8
1-70G	FG-27
1-77G	FG-10
1-332	F-690
1-360	FG-77
1-669	FG-28
1-690	FG-13
P-1	FG-49
P-12	F-324

Anti-Fatigue Mats

- Cushioned anti-fatigue mats reduces the risk of fatigue or static posture often caused by standing for long periods of time
- PharmaSystems® is proud to offer the following selection of anti-fatigue matting for use in dry environments
- Anti-fatigue mats are made to order. Additional lengths available with no shipping delays - contact our customer service team for options and pricing


How to select your anti-fatigue mats:

- Select by compound and thickness
- **Tabor Abrasion** – how well the surface withstands wear; the less weight lost the more durable the product (based on Federal Standard 191/1000g, CS-17 wheel; percentage lost refers to amount lost at 1000 cycles)
- **Coefficient of Friction** – degree of traction provided (OSHA recommends .5 in dry areas), based on ASTM F1677
- **Compression Deflection** - the softness or comfort level (ideal range is considered .20-.50); measured by percentage or inches @ 35 psi (Ψ)



Ribbed Anti-Fatigue Mats

- Economic comfort for light industrial applications
- Safety round edges
- Made-to-order lengths available, contact our customer service team
- Thickness: 3/8" (9 mm)
- Material: PVC Sponge
- Temperature Range: -17 - 43° C (0-110° F)
- Tabor Abrasion: 2% Lost
- Coefficient of Friction: Dry = .94 / Wet = .72
- Compression Deflection: ASTM D575/ .205 @ 35 Ψ



Item #	Size	Depth	Colour	Weight	UPC
5270	2 x 3'	3/8"	● Grey	4.45 lbs.	063636752705
5272	3 x 5'	3/8"	● Grey	10.1 lbs.	063636752729
5274	3 x 12'	3/8"	● Grey	26.14 lbs.	063636752743



Custom Mat Options!

Anti-fatigue mats to your specifications & imprinting options available!
Contact our customer service team to place your custom order today!

Endurable Mats


- Economic comfort for medium traffic areas
- Dual-layer construction outlasts lower quality blown sponges
- Safety beveled edges
- Made-to-order lengths available, contact our customer service team
- Thickness: 1/2" (13 mm)
- Material: Dual Density PVC Sponge
- Tabor Abrasion: 2% Lost
- Coefficient of Friction: Dry = .91
- Compression Deflection: ASTM D575/ .23 inches @ 35 Ψ

Item #	Size	Depth	Colour	Weight	UPC
45935	3 x 5'	1/2"	● Black	14.1 lbs.	063636745394
459312	3 x 12'	1/2"	● Black	40.1 lbs.	063636745936
45934	3 x 4'	1/2"	● Black	-	063636745967
459310	3 x 10'	1/2"	● Black	-	063636745905

Pharmacy Supplies

Traceable® Thermometers:

- Triple display simultaneously shows the current, minimum and maximum temperature reading
- Includes Traceable® certificate from ISO 17025 calibration laboratory, accredited by A2LA, valid for 1 year - recalibration is the customer's responsibility, contact us for assistance
- Alarm alerts when temperature rises or falls outside set point
- Low battery symbol indicates when batteries should be replaced
- Reads temperatures ranging from - 50° - 70° C (-58° - 158° F)
- Includes all necessary batteries and mounting equipment (flip stands, wall mounts, Velcro® and magnet strips)



Item #10369 is a dual bottle probe extension for item #10366

Traceable® Thermometers

- Reads the ambient room temperature as well as the external probe temperature using the in/out function
- Alarm programmable in 1° increments

Item #	Description	Probe Type	Dimensions (W x H x D)	UPC
14144	A Traceable® Refrigerator/Freezer Thermometer	5 mL Bottle Probe	2 3/4 x 4 1/4 x 3/4"	063636141448
10367	B Traceable® Refrigerator/Freezer Thermometer	Ambient Probe	2 3/4 x 4 1/4 x 3/4"	063636103675
10368	C Traceable® Refrigerator/Freezer Thermometer	10 mL Bottle Probe	2 3/4 x 4 1/4 x 3/4"	063636103682
17751	D Traceable® Jumbo Refrigerator/Freezer Thermometer	10 mL Bottle Probe	3 7/8 x 4 3/8 x 7/8"	063636177515
17752	Traceable® Jumbo Refrigerator/Freezer Thermometer	5 mL Bottle Probe	3 7/8 x 4 3/8 x 7/8"	063636177522
17750	E Traceable® Jumbo Refrigerator/Freezer Thermometer	Ambient Probe	3 7/8 x 4 3/8 x 7/8"	063636177508
14140	F Traceable® High Resolution Refrigerator/Freezer Plus™ Thermometer	10 mL Bottle Probe	3 1/4 x 2 1/5 x 5/8"	063636141400
10366	G Traceable® Dual Thermometer	2 10 ml Probes	3 1/2 x 1 3/4 x 1"	063636103668
10369	Traceable® Dual Bottle Probe Extension for #10366	2 Ambient Probes	3 1/2 x 1 3/4 x 1"	063636103668

New Pharmacy Recommendations

The following are our recommendations for new pharmacy opening orders.

Contact our customer service department for more information and details about your province's College of Pharmacy requirements, and for required reference books.

Starter Package | Item # ORD1

Item #	Item Description
2998	Filter paper
3002	Graduated glass cylinder 10ml
3006	Graduated glass cylinder 50ml
3008	Graduated glass cylinder 100ml
3016	Weighing paper
3017	Glass stirring rods 8" (3)
3091	Spatula rubber 4"
3042	Spatula stainless steel 4"
3028	Spatula stainless steel 6"
3029	Spatula stainless steel 8"
3030	Funnel (Plastic 55 mm) 2"
3031	Funnel (Plastic 100 mm) 4"
3073	Mortar & Pestle Glass Set 4 oz
3075	Mortar & Pestle Glass Set 8 oz
3076	Mortar Glass 16 oz
3086	Pestle Glass 16 oz
3088	Ointment Pad
3090	Weighing boats plastic 100/pkg
VIC-303	EJ Series Electronic Scale, 300 gram capacity
CAL200	200g Calibration Weight
17751	Traceable® Temperature Recorder - Jumbo Display

Deluxe Package | Item # ORD2 (includes items from ORD1)

Item #	Item Description
20050-12	Bottles Amber Glass, 50 ml 264/Case
200100	Bottles Amber Glass 100 ml 70/Case
200200	Bottles Amber Glass 250 ml 96/Case
200500	Bottles Amber Glass 500 ml 48/Case
2001000-01	Bottle Amber Glass 1L
20020	White Lined Caps, 100/Pkg (For #20050-12 & 200100)
20022	White Lined Caps, 150/Pkg (For #200200)
20024	White Lined Caps, 55/Pkg (For # 200500)
20026	White Lined Caps, 25/Pkg (For #2001000)
3014	Glass Graduated Cylinder, 250 ml
3015	Glass Graduated Cylinder, 500 ml
3371	Glass Graduated Cylinder, 1000 ml
3098	Glass Ointment Slab, 12 X 12", 1/4" Thick
7782	7 MI Sterile Dropper Bottles With Dropper Caps, 12/Case
7784	15 MI Sterile Dropper Bottles With Dropper Caps, 12/Case
D-50	Tablet Counting Tray And Spatula
BK-40	Warning Label System & 17 Labels
T1-1	Tape Rx 1" X 36 Yds., 1" Core
T1-3	Tape Rx 1" X 72 Yds., 3" Core
Td1	Tape Dispenser For 1" Wide Tape (1" And 3" Core)
PRP-3	Rx Receipts 100/Pad x 10 Pads
11194	Small Box For Rx Prescription Pads (Each)
10196	Large Box For Rx Prescription Pads (Each)
1976	Rx Prescription Pad File Folder (Double Fold), 100/PK
5272	Anti-Fatigue Mat Ribbed Grey, 3' x 5'

Item Number:	Qty:
<input type="checkbox"/> ORD1	<input type="checkbox"/>
<input type="checkbox"/> ORD2 (includes items from ORD1)	

Fax or Email Completed Order Form: 1.888.475.7155 | orders@pharmasystems.com

Pharmacy:

Address:

City:

Province:

Postal Code:

Date:

Name:

Phone:

Fax:

Email:

Wholesaler:

Account No.: