

Introducing Traceable® 7600 Smart Wi-Fi Connected Data Logger for Refrigerator and Freezer Temperature Monitoring

Features an interactive touch screen, smart calibrated sensor probes, open-door monitoring, remote alarm acknowledgement, and cloud connectivity through TraceableLIVE®.

Includes a NIST-traceable certificate of calibration from an A2LA-accredited ISO/IEC 17025:2017 laboratory.



Intuitive Touch Screen

Simple, one-time set up with touch-screen interface for quick access to monitoring data, settings, and alerts.



Smart Sensor Probes

New calibrated smart sensor probes. Keep the data logger, replace the sensors for more efficient recalibration.



Remote Monitoring

Acknowledge alarms remotely, and receive notifications for temperature excursions, open doors, and other critical events.



Enhanced Security

Secure access controls with optional PIN protection to prevent unauthorized changes.



Smart Sensor Calibration

Each sensor is calibrated instead of the data logger. When it's time to recalibrate, the data logger stays put, and you can swap out the sensor that's expired for one that is ready to go.

TraceableLIVE® Connectivity

Compatible with TraceableLIVE®, a cloud-based subscription software solution for connected data logging devices.

- ✓ Email, text, and push notifications
- ✓ Voice call escalations available with qualifying subscriptions
- ✓ Wi-Fi cloud-connected access to data
- ✓ Multiple users and multiple devices supported
- ✓ Meets 21 CFR Part 11 North American compliance and Annex 11 for Europe

Included Documentation

Traceable® products have an individually numbered Traceable Certificate that assures accuracy from an A2LA-accredited ISO/IEC 17025:2017 (1750.01) calibration laboratory.

How to Order

Available as complete kits, or purchase the monitor, probes, and sensors separately as needed.

Can support up to two of each - create a custom system for your individual needs.

Click each item number below to view specifications, pricing, and place your order online.



<u>10377</u>	Traceable® 7600 Smart Wi-Fi Data Logging Refrigerator/Freezer Thermometer with 1 bottle probe and 1 door sensor	
<u>10379</u>	Traceable® 7600 Smart Wi-Fi Data Logging Refrigerator/Freezer Thermometer with 1 bullet/ambient probe and 1 door sensor	
<u>10380</u>	Traceable® 7600 Smart Wi-Fi Data Logging Refrigerator/Freezer Thermometer with 1 vaccine probe and 1 door sensor	
<u>10371</u>	Traceable® 7600 Smart Wi-Fi Data Logger (only)	
<u>10372</u>	Smart Sensor Calibrated Temperature Bottle Probe (only)	
<u>10373</u>	Smart Sensor Calibrated Temperature Bullet/Ambient Probe (only)	
<u>10374</u>	Smart Sensor Calibrated Temperature Vaccine Probe (only)	
<u>10376</u>	Smart Sensor Calibrated Open-Door Sensor (only)	

Know What's Happening in Your Refrigerator or Freezer, No Matter Where You Are

TraceableLIVE® is a cloud-based subscription software solution for connected data logging devices that sends alerts to your Wi-Fi connected device when you need to take action.

- Alerts, emails, and push notifications sent to your smartphone or Apple Watch,* and voice call escalations* to your phone, for peace of mind that your samples are safe
- Wi-Fi cloud-connected, so your data is available whenever you need it
- Regular updates are automatically pushed to your system, keeping it maintenance free, up to date, and secure
- Quick to set up and easy to use, keeping you focused on what matters most
- Add multiple users and multiple devices for robust lab monitoring



TraceableLIVE® Activation and Subscription Plans

TraceableLIVE® service comes in three levels:

	Free Account	Standard Account	Premium Account
Pricing	Free	1-Time Activation Fee, plus Subscription (one device, per month)	1-Time Activation Fee, plus Subscription (one device, per month)
Users	1 (only 1 admin)	Unlimited admins and users	Unlimited admins and users
Max Devices	3	No limit (each device must have subscription)	No limit (each device must have subscription)
Data Storage	Rolling 30 days	Unlimited	Unlimited
Email Notifications	✓	✓	✓
Push-to-App Notifications	✓	✓	✓
Text Notifications	–	–	✓
Escalation Voice Call Notifications	–	–	✓

Purchase Subscription Activation Cards

Click each item number to view specifications, pricing, and order online.

1 – Year Standard
Part Code 10360

2 – Year Standard
Part Code 10361

1 – Year Premium
Part Code 10362

2 – Year Premium
Part Code 10363