

Are you Ready for Vaccination Season?

Ordering Vaccine Supplies is now easier than ever with PharmaSystems:

- ✓ In-Demand Supplies
- ✓ Industry-Leading Service
- ✓ Same-Day Order Processing



Safety Needles

- Single use, standard hypodermic needle with a protective arm that is used to enclose the needle after use to protect the healthcare professionals from needlestick injury
- Patented locking mechanism ensures needle is secure and completely locked in protective arm
- Compatible with Luer Lock syringes at right for all vaccinations, including Flu and Covid products



Product Code	Description	Box Qty
V0205*	25G x 1" Adult Needle	100
V0206*	25G x 1-1/2" Higher BMI Population	100
V0207*	25G x 5/8" Adolescents/Infant	100

[VIEW ONLINE](#)

Luer Lock Syringes

- 1ml Low Dead Space syringes minimize "dead space," the area that retains medication even after it has been injected - this improves the efficiency of medication delivery by reducing the amount of medication wasted with every administration
- High-transparency barrel
- Silicone lubrication ensures smooth plunger movement
- Safe plunger backstop prevents accidental withdrawal of plunger

Product Code	Description	Box Qty
V0201*	1mL, Luer Lock Syringe	100
V0203	3mL, Luer Lock Syringe	100



[VIEW ONLINE](#)

*Low Dead-Space

Consolidate your Vaccine Clinic Supplies through a Single Source

Connect with our bilingual customer service team:

Order Online: www.pharmasystems.com

Call: 905.475.2500 (Toll free: 1.888.475.2500)

Fax: 905.475.7155 (Toll free: 1.888.475.7155)


Email: orders@pharmasystems.com



Register your Online Account:

- ✓ **Easy Ordering** Quick & efficient order & re-ordering
- ✓ **Purchasing Efficiencies** Consolidate supply orders
- ✓ **Same-Day Order Processing**
- ✓ **Category Expertise** Connect with our industry-leading customer service team
- ✓ **Select Preferred Billing Method** Select to dropship through your wholesaler, or direct-to-store order

View Vaccine Products:



[REGISTER](#)

Quick Order Form



Vaccine Essentials



Please Email or Fax Completed Form to: orders@pharmasystems.com | Fax: 1.888.475.7155

Part Code	Description	Pack Size	Enter Order Quantity
Safety Needles			
<u>V0205*</u>	25G x 1" Adult Safety Needles	100	
<u>V0206*</u>	25G x 1-1/2" Higher BMI Population Safety Needles	100	
<u>V0207*</u>	25G x 5/8" Adolescents/Infant Safety Needles	100	
Luer Lock Syringes			
<u>V0201*</u>	1ml Luer Lock Syringe LDS (low dead space), without Needle	100	
<u>V0203</u>	3ml Luer Lock Syringe, without Needle	100	
Sharps Disposal			
<u>9604-01</u>	Sharps Container, 1 Quart	1	
<u>9601-01</u>	Sharps Container, 5 Quart (horizontal drop)	1	
<u>9601-06</u>	Sharps Container, 5 Quart (horizontal drop)	20	
Gloves			
<u>GNSMS</u>	Nitrile Gloves (Blue), Small	100	
<u>GNSMM</u>	Nitrile Gloves (Blue), Medium	100	
<u>GNSML</u>	Nitrile Gloves (Blue), Large	100	

*Low Dead-Space. Learn more here: pharmasystems.com/vaccine-clinic-supplies

PHARMACY NAME:		CONTACT NAME:	
ADDRESS:		PHONE:	
CITY:		FAX:	
PROVINCE:	POSTAL CODE:	EMAIL:	