

Supporting **Every Step** of Veterinary Care

Products, tools & solutions for
everyday animal care.





PharmaSystems is proud to support the day-to-day needs of veterinary practices across Canada.

This catalogue features a curated selection of essential clinic and pharmacy supplies designed to enhance efficiency, organization, and patient care.

Discover high-demand veterinary clinic supplies, dispensing and packaging solutions, examination room essentials, prescription labeling products, safety and infection prevention supplies, and operational necessities trusted by veterinary hospitals and animal clinics nationwide.

Canadian veterinary professionals can rely on these quality products to help streamline workflows, support preventive care, and deliver exceptional patient outcomes in today's busy veterinary environment.


For our full selection of veterinary care products and pharmacy supplies, log in to our online portal today **here**.

.....

Please Note: We have developed this curated catalogue to showcase our most requested veterinary supplies. Please refer to our website for our entire product catalogue and for the most up-to-date product information. Product images are for illustrative purposes only. Please refer to the product description and specifications for complete details.

Terms: Net 30 days or as per your wholesaler

Returns: All returns to PharmaSystems must be authorized.



Create an account to place your order today!

Benefits

- View pricing: Simply register your account
- Place orders online through our online shopping cart, or by downloading a blank order form (return by fax or email)
- Be the first to know about new products
- Save on deals & promotions: Learn about new product promotions and take advantage of deals & clearance sales
- Stay informed: Keep up to date with changing Canadian Pharmacy Regulations and learn about product recommendations.

How to order:

Phone in your order: 905.475.2500 | Toll free: 1.888.475.2500

Fax in your order: 905.475.7155 | Toll free: 1.888.475.7155

Email your order: orders@pharmasystems.com

Use the order form found on the last page

Our bilingual customer service team is readily available to help you with any of your pharmacy supply requests!

T: 905.475.2500 | Toll free: 1.888.475.2500

F: 905.475.7155 | Toll free: 1.888.475.7155

Address: 151 Telson Road Markham Ontario Canada L3R 1E7

Email: orders@pharmasystems.com

Create an account here: www.pharmasystems.com



TABLE OF CONTENTS

05 BAGS

- Reclosable Bags
- Prescription Bags
- IV Covers & Bags

10 BLOOD COLLECTION

13 ORGANIZATION & STORAGE

- Boxes & Bins
- Caddies
- Rx Fill Trays

21 COMPOUNDING & DISPENSING

- Beakers & Cylinders
- Mortar & Pestles
- Spatulas
- Bottles & Vials

32 CRUSHERS

35 MEDICATION DELIVERY

- ENFit Supplies
- Oral Dispensers
- Bottles & Adapters
- Medication Cups
- Syringes & Needles

47 PROTECT


- PPE
- Infection Control
- Spill Kits
- Sharps Containers

54 MONITOR

- Thermometers
- Cold Storage

57 IDENTIFY & SECURE

- Tape
- Tamper-Evident Seals
- Labels



BAGS

Reclosable Bags
Prescription Bags
IV Covers & Bags



Resealable Bags

[Order Online](#)
[View More Products](#)

Features

- Thickness: 3-mil
- Single-track zipper closure allows quick access while helping keep contents protected between uses.

Item No.	Size	Colour	Qty
7567	2.5 x 9"	Amber	100
7567M	2.5 x 9"	Amber	1000
7561	2 x 3"	Amber	100
7570	3 x 5"	Amber	100
7570M	3 x 5"	Amber	1000
7573	4 x 6"	Amber	100
7573M	4 x 6"	Amber	1000
7575	6 x 8"	Amber	100
7575M	6 x 8"	Amber	1000
7576	8 x 8"	Amber	100
7576M	8 x 8"	Amber	1000
7588	8 x 14"	Amber	100
7588M	8 x 14"	Amber	1000
7578M	12 x 12"	Amber	100
7578M	12 x 12"	Amber	1000
M4G3X5	3 x 5"	Clear	1000
MG2X3	3 x 3"	Clear	1000
MG4X4	4 x 4"	Clear	1000
MG4X6	4 x 6"	Clear	1000

TIP

Use amber bags for medications that require protection from light to help maintain product stability during storage and transport.



UV Zip-It Bags, Amber

[Order Online](#)[View More Products](#)

Designed to protect light-sensitive medications, UV Zip-It® Bags provide secure storage with convenient visibility and easy access.

Features


- Easy to access and provides sturdy storage
- Latex-free
- Protects light sensitive medications - amber colour bags absorb 92% to 98% of light in the UV range - the dark brown and black colour bags absorb 99.9% of light in the UV range
- Lighter colours allow the nurse to read a label or check the level of an IV fluid through the UVLI bag

Item No.	Size	Qty
AUZ35	3 x 5"	1000
AUZ39	3 x 9"	1000
AUZ46	4 x 6"	1000
AUZ58	5 x 8"	1000
AUZ69	6 x 9"	1000
AUZ814	8 x 14"	1000
AUZ912	9 x 12"	1000

TIP



The slit-top design allows fast access to contents while maintaining a clean, professional appearance.



Prescription Bags

[Order Online](#)
[View More Products](#)

Hanging Prescription Bags

- Simplify prescription organization and retrieval
- Made of tough, clear 4 mil plastic for extreme durability
- Fits on any rod with a diameter of 1 cm or smaller

Wall Rack

- Simple solution to prescription organization
- Saves time, improves accuracy and enhances customer satisfaction
- Aluminum frame is 3 ft. wide and holds five 1 cm thick rods
- Each rod holds approximately 8-10 bags per foot or as many as 150 bags per rack
- Easy assembly, all hardware included
- Measures: 36" W x 57.5" H x 7.5" D
- Use with hanging bags and ID tags (sold separately)



Item No.	Description	Size	Qty
17540	Hanging Prescription Bag	10.5 x 7.5"	10
17541	Hanging Prescription Bag	12 x 9"	10
17542	Hanging Prescription Bag	14 x 12".5"	10
17543	Wall Rack	36" W x 57.5" H x 7.5" D	1 frame, 5 rods
17544	Green ID Tags	-	25

TIP

Organize hanging bags by patient name, prescription status, or pickup time to streamline operations and enhance customer service.



IV Covers & Bags

[Order Online](#)
[View More Products](#)

Features

- Bags are transparent, but will mask many differences in IV fluid color.
- Thickness: 1.5-mil - 1.8 mil

Item No.	Size	Description	Tint	Qty
7589	8 x 14"	IV Cover, slit at sealed end	Light	100
7587	6 x 10"	IV Cover, slit at sealed end	Amber	100
7587D	6 x 10"	IV Cover	Dark	100
7589D	8 x 14"	IV Cover, slit at sealed end	Dark	100
7590D	8 x 14"	IV Cover	Dark	100
19643	12 x 18"	IV Cover	Amber	100
7583	8 x 14"	IV Cover, slit at sealed end	Green	100
7584	6 x 10"	IV Cover, slit at sealed end	Green	100
7590L	8 x 14"	IV Cover	Amber	100
7579	5 x 8.5"	IV Cover, slit at sealed end	Green	100
590-UT300	2.5" x 300'	IV Tubing Cover, opened end	Amber	1 Roll
7617	2.5 x 72"	UV Light Covers	Amber	10 Strips


DID YOU KNOW?

Dark Amber covers block 95% of UV light, while Light/Amber and Green covers block 90% of UV light. They are perfect for dispensing IV admixtures that contain light-sensitive drugs.



BLOOD COLLECTION

Test Tube Storage
Phlebotomy Trays



Test Tube/Syringe Rack | 7405

[Order Online](#)

[View More Products](#)

Features

- Metal construction has beige enamel finish.
- Holds 32 syringes, test tubes or blood tubes.
- Each of the 32 sections measures 5/8" x 5/8"
- Overall Dimensions: 6"W x 4-5/8"H x 3"D



TIP

Assign specific rows for different syringe sizes or sample types to streamline workflow and reduce retrieval time.



Specimen Caddy | 5607

[Order Online](#)


[View More Products](#)

Features

- Maximize Storage - 6 compartment bottom drawer that keeps supplies safe and out of patients' and visitors' view.
- Quick Access - This Specimen Caddy offers quick and easy access to blood collection supplies.
- Safety Latch - Drawer features an easy to open safety latch that eliminates accidental opening.
- Includes: One sharps bracket that accommodates a sharps container and 2 test tube/syringe racks that can hold 32 tubes per rack.
- Dimensions:
 - Outside - 14^{3/4}"W x 6"H x 9-13/16"D
 - Height Including Handle - 6-13/16"

TIP

Designate each compartment for a specific supply type to streamline specimen collection and make restocking faster and more efficient.



Phlebotomy Tray

[Order Online](#)
[View More Products](#)

Features

- Versatile - Features a set of two economical, disposable, compartmentalized inserts specifically designed for blood collection.
- Sort Contents - Easy to organize supplies such as prep pads, needle holders, lancets, etc. for quick retrieval.
- Easy to Clean - Wash as needed and replace the insert periodically for a fresh start.
- Carrying Handle - Ratcheting handle with soft rubber grip rotates out of the way for easy access.
- Includes: Sturdy plastic tray comes with its own 16mm test tube rack and two disposable Phlebotomy Tray Inserts (17466). Supplies are not included with tray.
- Dimensions:
 - Outside - 16-7/8"W x 5-7/8"H x 9¼"D
 - Height Including Handle - 11¼"

Inserts:

- Efficient - These disposable tray inserts improve patient care by providing a more consistent, efficient, and convenient method of maintaining clean Phlebotomy Trays (sold separately).
- Maximize Space - Tray inserts include eight compartments per insert. Use 16mm Tube Racks with this insert to organize tubes.
- Use with Phlebotomy Tray (17465).

Item No.	Description	Qty
17465	Phlebotomy Tray with 2 Inserts	1
17466	Extra Tray Inserts	4

TIP

Assign each compartment a specific purpose to create a standardized setup that speeds up supply retrieval and simplifies restocking.



ORGANIZATION & STORAGE


Boxes & Bins

Caddies

Rx Fill Trays

Tilter Bins

Totes



Organization Bins

[Order Online](#)
[View More Products](#)

Divider Boxes

- Smart Design: Recessed stacking ledges protect contents of container below
- Durable Plastic: Resistant to most solvents, chemicals, stress cracks
- Autoclavable: 250°F (121°C) for two minutes.
- Optional Accessories - Dividers, lids and label holders are sold separately
- Customizable: Plastic dividers dividable every 1¼"
- Color Options: Blue, Gray, Red or Light Blue. Specify color when ordering.

Lids


- Lids include lift-holes providing an area to secure the box with security seals (sold separately)
- Semi-clear snap-on lids fit over the top edges of the box
- Clear deluxe lids fit on inside rim of boxes
- Snap on Lids: 1701, 1706, 1716
- Deluxe Lids: 1730, 1730-01, 1731, 1731-01, 1732

Item No.	Outside Dimensions (L X W X H)	Colour	Compatible Lids
1700	10.875 x 8.25 x 5"	● ● ●	1701, 1730
1705	16.5 x 10.875 x 6"	● ● ●	1706, 1731
1710	16.5 x 10.875 x 8"	● ● ●	1706, 1731
1715	22.375 x 17.375 x 8"	● ● ●	1716, 1732
1723	16.5 x 10.875 x 3.5"	● ● ●	1731-01
1724	10.875 x 8.25 x 2.5"	● ● ●	1730-01
1725	10.875 x 8.25 x 3.5"	● ● ●	1730-01
1891	16.5 x 10.875 x 4"	● ● ●	1706, 1731
1892	22.375 x 17.375 x 10"	● ● ●	1716, 1732
1893	22.375 x 17.375 x 4"	● ● ●	1716, 1732
1894	22.375 x 17.375 x 6"	● ● ●	1716, 1732

TIP

Combine adjustable dividers, labels, and color-coded bins to create a highly organized storage system that makes inventory identification and retrieval faster and more efficient.

View dividers [here](#).



Tilters® Tilt Front Bins

[Order Online](#)
[View More Products](#)

Features

- Modular bins: Stack and customize to fit your storage and organization requirements
- Front opens at a 45 degree angle for easy access
- Clear fronts provide easy visual access to inventory
- Available in beige or white. Specify colour when ordering.

Item No.	Overall Dimensions (L X W X H)	Inside Dimensions (L X W X H)	Compartments	Colour
1203A	23.625 x 4.5 x 3.625"	3.125 x 2.75 x 2.625"	6	○○
1205E	23.625 x 6.5 x 5.25"	4 x 4.5 x 3.75"	5	○○
1207	23.625 x 8.125 x 6.62"	5 x 5.25 x 5"	4	○○
1210A	23.625 x 9.5 x 7.75"	7.75 x 23.625 x 9.5"	3	○○
1215	23.625 x 13.75 x 9.62"	11.875 x 23.625 x 13.875"	2	○○
1213	11.812 x 13.87 x 12.25"	10.5 x 9.375 x 8.25"	1	○○



TIPS

The 45° tilt-open design provides quick access to contents while helping keep items securely contained.

Modular storage bins can help maximize vertical storage space, making them ideal for crowded work areas.



STORE | BINS

Rx Fill Tray

[Order Online](#)

[View More Products](#)

Features

- Better Rx Management: Keep individual prescription components together during the entire filling process.
- Space Saving: Trays with easy-grip handles nest when empty.
- Built-in Holder: Tray cutout allows use with a standard spring clip to hold important documentation.
- Dimensions: 6-7/8"W x 2½"H x 10"D
- Sold individually

Item No.	Colour
17487B	Blue
17487C	Clear
17487R	Red
17487Y	Yellow
17487G	Green
17487W	White



TIP

Nest empty Rx Fill Trays to maximize cabinet space and keep prescription workflows organized.



Heavy Duty Storage Bins

[Order Online](#)[View More Products](#)

Features

- Shield Supplies: Transparent lids keep contents dust-free and protected.
- Easy Content ID: Label slots hold fanfold, roll or laser labels
- Sturdy: Reinforced sides keep bins from spreading.
- Compatible: Waterproof, oil and chemical resistant.
- Weight Capacity: 30 lbs.
- Dimensions:
Outside: 4 $\frac{1}{8}$ "W x 4"H x 10 $\frac{7}{8}$ "D
Inside: 3-13/16"W x 3 $\frac{3}{4}$ "H x 10 $\frac{1}{4}$ "D
Base Length = 9"

Item No.	Colour
19322K	Black
19322B	Blue
19322PP	Purple
19322R	Red
19322Y	Yellow
19322G	Green
19322W	White



TIP

Use labeled storage bins with transparent lids to quickly identify and protect inventory while keeping supplies organized, dust-free, and within easy reach.



Cardboard Shelf Caddies

[Order Online](#)[View More Products](#)

Features

- Budget Friendly: Our most economical bins ship flat for reduced shipping costs.
- Convenient: Write-on surface for quick labeling.
- Disposable
- Compact Storage: Ships flat

Item No.	Dimensions (W x H x D)	Qty
7475	12 x 4.5 x 12"	50
7458	2 x 4.5 x 12"	50
7460	4 x 4.5 x 12"	50
7465	6 x 4.5 x 12"	50
7470	9 x 4.5 x 12"	50

TIP

Cardboard bins can be recycled. When a bin is dirty, simply replace it. No need to clean!



STORE | BINS

Smart Store Tote with Lid

Item # 13204C

Order Online

View More Products

Features

- Space Saving: Totes will nest when empty for space saving storage.
- Durable: Withstands temperature ranges from -40°F to 248°F (-40°C to 120°C)
- Material: Polypropylene
- Construction:
 - Double latches to ensure that lids remain with totes
 - Sides feature built-in carrying handles for easy carrying
 - Reinforced heavy-load bottoms
- Capacity: 15-Liter
- Dimensions: 15¾"W x 7¼"H x 11-13/16"D



TIP

Use durable totes to safely transport and organize inventory, then nest them when empty to save valuable storage space.



Shelf Bins

[Order Online](#)[View More Products](#)


Features

- Sort Supplies: Plastic bin fits neatly inside cart drawers and in tight spaces to organize smaller quantities of medication and supplies.
- Improve Organization: Add up to nine Small Bin Cups or four Large Bin Cups (sold separately) to help sort and organize contents.
- Easily ID Contents: Insert labels in the bin label holder for quick supply identification. Use with the Label Extension Holder to improve visibility. (Labels and Label Extension Holder sold separately)
- Bins nest when not in use.
- Dimensions:
Outside: 6-5/8"W x 4"H x 11-5/8"D
Inside: 5½"W x 4"H x 10¼"D

Item No.	Colour	Qty
1445R	Red	12
1445B	Blue	12
1445C	Clear	12
1445Y	Yellow	12

TIP

Use shelf bins with labels and bin cups to keep medications and supplies neatly sorted, easy to identify, and readily accessible in drawers, carts, and storage areas.



COMPOUNDING & DISPENSING SUPPLIES

Beakers

Bottles & Vials


Cylinders

Mortar & Pestles

Pill Counting Trays

Spatulas

Weighing Boats



Glass Beakers

[Order Online](#)
[View More Products](#)

Features

- Scale: White enamel
- Clean Results: Spout for easy pouring
- Compatible: Safe for use on hot plates

Item No.	Volume	Increments	Qty
3372	30 mL	10 mL	1
3024	50 mL	10 mL	1
3025	100 mL	10 mL	1
3373	150 mL	25 mL	1
3026	250 mL	25 mL	1
3374	400 mL	25 mL	1
3027	600 mL	50 mL	1
3375	1000 mL	50 mL	1

TIP

Pour slowly with the spout pointed directly at the receiving container, and keep the beaker close to the opening to reduce splashing and drips.



Glass Graduated Cylinders

[Order Online](#)
[View More Products](#)


Features

- Made from high-quality borosilicate glass
- Single metric scale
- Has an easy-pour spout with large circular base
- Rated at 20°C
- Not recommended for autoclaving; for autoclaving please refer to 1634AA series [here](#)

Item No.	Volume	Increments	Qty
3002	10 mL	0.2 mL	1
3009	25 mL	0.5 mL	1
3006	50 mL	1 mL	1
3008	100 mL	1 mL	1
3014	250 mL	2 mL	1
3015	500 mL	5 mL	1
3371	1000 mL	10 mL	1

TIP

Select the smallest graduated cylinder that accommodates your required volume for greater measurement precision.



Flint Glass Mortar & Pestles

[Order Online](#)
[View More Products](#)

Features

- Convenient: Use for preparing medications for patients or compounding recipes
- Easy to Clean: Can be sterilized, cleaned with alcohol or washed with soap and water
- Note: Not recommended for use with drastic temperature changes
- Flint glass models are ideal for compounding liquids, ointments & pastes
- Quantity: Sold individually

Item No.	Description	Size
3070	Flint Glass Mortar	2 oz. (60 mL)
3080	Flint Glass Pestle	2 oz. (60 mL)
3071	Flint Glass Mortar & Pestle Set	2 oz. (60 mL)
3072	Flint Glass Mortar	4 oz. (120 mL)
3082	Flint Glass Pestle	4 oz. (120 mL)
3073	Flint Glass Mortar & Pestle Set	4 oz. (120 mL)
3074	Flint Glass Mortar	8 oz. (240 mL)
3084	Flint Glass Pestle	8 oz. (240 mL)
3075	Flint Glass Mortar & Pestle Set	8 oz. (240 mL)
3076	Flint Glass Mortar	16 oz. (480 mL)
3086	Flint Glass Pestle	16 oz. (480 mL)
3382	Flint Glass Mortar & Pestle Set	16 oz. (480 mL)

TIP

Allow glass equipment to reach room temperature before cleaning or sterilizing to help prevent damage from sudden temperature changes.



Compounding Accessories

[Order Online](#)
[View More Products](#)

Ointment Pad

- 12" square
- 100 sheets per pad

Stainless Steel Spatula

- Versatile - Flexible stainless steel blade with comfortable wooden handle
- Easy Cleaning - Use soap and water



Plastic Weighing Boats

- Inexpensive containers for weighing liquids or solids.
- Resilient polystyrene construction
- Non-stick surface, rugged enough for mixing
- Stack and nest easily

Item No.	Description	Dimensions	Qty
3088	Ointment Pad	12" x 12"	100/pad
3042	Spatula Stainless Steel	4" Blade	1
3043	Spatula Stainless Steel	5" Blade	1
3028	Spatula Stainless Steel	6" Blade	1
3029	Spatula Stainless Steel	8" Blade	1
3429	Spatula Stainless Steel	10" Blade	1
3430	Spatula Stainless Steel	12" Blade	1
3420	Plastic Weighing Boat Micro	4.4 x 4.4"	100
40100	Plastic weighing Boat, Small	4.5 x 4.5"	100
3090	Plastic weighing Boat, Medium	8.3 x 8.3"	100
3092	Plastic weighing Boat, Large	14 x 14"	100

TIP

Flexible stainless steel spatulas help achieve smoother blending and more uniform incorporation of ingredients during compounding.



Plastic Dropper Bottles, Sterile

[Order Online](#)
[View More Products](#)

Engineered for sterile compounding applications, Plastic Dropper Bottles combine secure storage with accurate drop-by-drop dispensing.

Features

- **Sterility Assurance:** Ideal for the preparation of sterile ophthalmic solutions. Bottle exterior is sterile until the package is opened; the inside is sterile until the cap is removed.
- **All-in-One System:** Includes snap-in dropper tips and screw on caps.
- Radiation sterilized
- USP 797 Compliant

Item no.	Size	Colour	Amount released per drop	Qty
7780	3 mL	Clear	0.0238 ml/g	12
7780-02	3 mL	Clear	0.0238 ml/g	1
7782	7 mL	Clear	0.035 ml/g	12
7782-01	7 mL	Clear	0.035 ml/g	1
7784	15 mL	Clear	0.035 ml/g	12
7784-01	15 mL	Clear	0.035 ml/g	1
7786	30 mL	Clear	0.035 ml/g	12
7786-01	30 mL	Clear	0.035 ml/g	1
10079*	30 mL	Amber	-	12



TIP

Amber bottles help protect light-sensitive contents from UV and visible light exposure.

Sterile Vials

[Order Online](#)
[View More Products](#)

- Type 1 Glass USP Borosilicate
- Latex Free Chlorobutyl Stopper
- Aseptically Assembled in an ISO 5 (class 100) Cleanroom
- USP Tested for Sterility and Endotoxins
- Amber glass vials are specifically designed to protect the sterility and integrity of light-sensitive substances

Item No.	Capacity	Colour	Pack Qty	Case Qty
20837	2 mL	Clear	25	-
7761-03	5 mL	Clear	25	100
7762-03	10 mL	Clear	25	100
7766-03	20 mL	Clear	16	-
7763-03	30 mL	Clear	25	100
7765-03	50 mL	Clear	25	-
7764-06	100 mL	Clear	25	-
19806	2 mL	Amber	25	-
19813	5 mL	Amber	25	-
18702	10 mL	Amber	25	-
18703	30 mL	Amber	25	-
18704	50 mL	Amber	25	-
19817	100 mL	Amber	25	-



DID YOU KNOW?

Amber sterile vials are specifically designed to help protect light-sensitive medications while maintaining the sterility and integrity of the contents.



Glass Dropper Bottles

[Order Online](#)
[View More Products](#)

Features

- Flexible: Bottles provide usage flexibility with optional accessories such as glass droppers (19456-19461) and tamper evident caps (10067, 10068).
- Glass Droppers: Ideal for dispensing topical treatments when used together with Amber Glass Bottles

Item No.	Description	Size	Qty	Compatible With:
19450	Amber Glass Dropper Bottles	5 ml	340	19456
19451	Amber Glass Dropper Bottles	10 ml	255	19457
19452	Amber Glass Dropper Bottles	20 ml	216	19458
19453	Amber Glass Dropper Bottles	30 ml	168	19459
19454	Amber Glass Dropper Bottles	50 ml	99	19460
19455	Amber Glass Dropper Bottles	100 ml	68	19461
19456	Glass Droppers	-	100	19450
19457	Glass Droppers	-	100	19451
19458	Glass Droppers	-	100	19452
19459	Glass Droppers	-	100	19453
19460	Glass Droppers	-	100	19454
19461	Glass Droppers	-	100	19455
10067	Tamper Evident Caps, White	-	100	All the above bottles
10068	Tamper Evident Caps, Black	-	100	All the above bottles

TIP

Amber Type 3 pharmaceutical glass bottle is ideal for dispensing topical treatments and other UV light-sensitive liquid medications.



Plastic Oval Medicine Bottles

[Order Online](#)

[View More Products](#)

Features

- USP Compliant: Child-resistant cap conforms to USP "Tight" standards.
- Better Accuracy: Easy-to-read graduations
- Label Ready: Flat back panel for easy label placement
- Colour: Amber
- Includes caps

Item No.	Size	Bottle Opening	Qty
19489	2 oz.	24 mm	100
19490	4 oz.	24 mm	100
19491	8 oz.	24 mm	100
20833	16 oz.	28 mm	50



TIP




Amber bottles help protect light-sensitive liquid medication from UV and visible light exposure.

Tablet Counters

[Order Online](#)
[View More Products](#)

Fretwell Triangle

- Precision-made scientific instrument
- Count regular shaped tablets with absolute accuracy and speed
- Simply count the number of rows to determine tablet count - unique legend determines count
- Tablets glide easily on the stainless steel triangle
- Easy to clean



Pill Counters

- Versatile - High quality counters designed specifically for right and left-handed users!
- Essential Tool - Pharmacists and pharmacy technicians can count and dispense pills from stock bottles to dispensing vials without touching hands.
- Easy Handling - Includes counting spatula with ergonomic grip.
- Stay Connected - Tray features a docking slot that stores the spatula with the tray.
- Hinged Lid Covering - The dispensing chute latches when closed to prevent accidental opening while filling vials.



Item No.	Description	Qty
FT1	Fretwell Triangle	1
D-50	Pill Counter With Spatula	1
7147	Tablet Counter, Right-Handed	1



DID YOU KNOW?

The Fretwell Triangle's unique row-counting design enables fast, accurate tablet counting, helping streamline dispensing tasks in busy pharmacy environments.



PREPARE | COMPOUNDING & DISPENSING


Medication Envelopes | 8020

Order Online

View More Products

Features

- Cost Effective - Sealable envelope for small quantities of medication.
- Write-On Area - For patient information and instructions.
- Quantity - 250 per package
- Dimensions - 3½"W x 2¼"H



DID YOU KNOW?

Complete all labeling before filling the envelope to help reduce dispensing errors and improve workflow efficiency.

CRUSHERS

Pill Splitters

Silent Knight

Pill Crusher Bags



CRUSHERS

Pill Splitters & Crushers

[Order Online](#)[View More Products](#)

Pill Splitter

- Splitter is designed to hold pills securely in place for precise splitting
- Easy to use, safe and convenient
- Blade is made out of surgical stainless steel
- Built-in storage container is designed to catch pills

Pill Crusher

- Turns tablets into powder effortlessly
- Dual storage compartments
- Cap is threaded for extra security

Item No.	Description	Qty
17043	Pill & Tablet Splitter	1
7474	Deluxe Pill & Tablet Splitter	1
7338	Pill Crusher	1



TIP

Always confirm that a medication is suitable for splitting or crushing before altering its form.



Silent Knight® & Crusher Bags

[Order Online](#)
[View More Products](#)

Crusher Pouches

- Secure Crushing - Tablets stay contained during crushing. (Do not crush enteric-coated tablets.)
- Convenient Storage - Compatible with Silent Knight® built-in pouch storage.
- Durable Design - Holds up to 5 tablets; helps reduce cross-contamination
- Thickness - 6.25 ± 0.5 mil
- Dimensions - 2"W × 4⅞"H

Silent Knight Tablet Crusher

- Powerful Performance - Quietly crushes tablets into a fine powder.
- Easy Operation - Requires minimal effort.
- Contained Crushing - Uses pouches to hold tablets during crushing.
- Built-In Storage - Holds crusher pouches (7414-31).
- Easy to Clean - Compatible with alcohol and most disinfectants.
- Dimensions - 10½"W × 4½"H × 3"D



Item No.	Description	Qty
7414-10	Tablet Crusher Pouches	100
7414-31	Tablet Crusher Pouches	1,000
7424-01	Silent Knight® Pill Crusher	Includes 50 pouches

DID YOU KNOW?

Crushing tablets directly within a Silent Knight® pouch helps contain medication, reduce product loss, and simplify cleanup



MEDICATION DELIVERY

Enteral Feeding Supplies

Oral Dispensers

Ampule Breakers

Dispensing Pins

Bottles & Adapters

Syringes

Needles

Luer Lock Connectors



ENFit® Supplies

[Order Online](#)
[View More Products](#)

ENFit Luer Lock Connectors

- Female/Male Compatible - Provides a connection between ENFit feeding tubes and luer lock syringes
- Simplifies the process of flushing ENFit feeding tubes with sterile saline.

ENFit Oral Stepped Bottle Adapter

- Stepped bottle adapter with male ENFit® connector and Oral connector
- For use with bottle sizes with an inner diameter of 16-28mm
- Intended to use for the lifetime of the bottle – once the bottle is empty, discard the adapter

Clear ENFit Syringes

- Patient Safety - ENFit syringes support ISO 80369-3 compliance.
- Prevents Misconnections - Twist-to-connect tip fits enteral devices only.
- Versatile - Suitable for feeding, flushing, hydration, aspiration, and medication delivery.
- Easy to Read - Clear 1 mL graduations improve dosing accuracy.
- Transition Ready - Includes adapter for use with legacy feeding tube sets
- Single Patient Use - One syringe per patient, per dispensing.

Item No.	Description	Capacity	Qty
20387	ENFit Luer Lock Connectors	-	100
21052	ENFit Oral Stepped Bottle Adapter	-	1
19800-500	Clear ENFit Syringes	5 ml	500
19801-500	Clear ENFit Syringes	10 ml	500

TIP

Not sure which bottle adapter will fit? The ENFit Oral Stepped Bottle Adapter is universal and fits a wide range of bottle sizes.



KwikDose™ Oral Dispenser

[Order Online](#)
[View More Products](#)

Designed for neonatal, pediatric, and adult use, KwikDose™ Oral Dispensers deliver safe and accurate oral dosing, while preventing wrong-route administration.

Features

- Self-righting caps included (will not fit luer needle)
- Unique tip design prevents wrong-route administration
- Chemically inert to provide stability for drug storage and transport
- Enlarged numbers for quick and accurate volume identification
- Includes distinct bilingual labeling “For Oral Use Only”



Amber Dispensers

Item No.	Capacity	Increments	Qty
010200	0.5 mL	0.01 mL	100
0200	0.5 mL	0.01 mL	500
010201	1 mL	0.01 mL	100
0201	1 mL	0.01 mL	500
010203	3 mL	0.01 mL	100
0203	3 mL	0.01 mL	500
010205	5 mL	0.2 mL	100
0205	5 mL	0.2 mL	500
010210	10 mL	0.2 mL	100
0210	10 mL	0.2 mL	500
010220	20 mL	1 mL	100
01221	60 mL	1 mL	75

Clear Dispensers

Item No.	Capacity	Increments	Qty
011200	0.5 mL	0.01 mL	100
01200	0.5 mL	0.01 mL	500
011201	1 mL	0.01 mL	100
01201	1 mL	0.01 mL	500
011203	3 mL	0.01 mL	100
01203	3 mL	0.01 mL	500
011205	5 mL	0.2 mL	100
01205	5 mL	0.2 mL	500
011210	10 mL	0.2 mL	100
01210	10 mL	0.2 mL	500
011220	20 mL	1 mL	100
01222	60 mL	1 mL	25



TIP

Amber oral dispensers helps protect light-sensitive medications from UV and visible light exposure.

DELIVER | MEDICATION DELIVERY

DoseLine Labels

Item No. 42910

.....

[Order Online](#)

[View More Products](#)

Exclusively available through PharmaSystems!

Features

- Clear sticker with black key line
- Designed for use on oral liquid dispensers, or other graduated containers
- Made of clear poly film, latex free
- Labels measure 1-9/16" x 3/8" (40 x 9.5 mm)
- Conveniently packaged in rolls of 500, 2 rolls per box
- Packaged in a convenient dispensing blister, compatible for use with PharmaSystems' label dispensers



TIP

Apply DoseLine Labels to oral syringes or graduated containers to create a clear, consistent dosing reference that helps improve measurement accuracy and makes liquid medication administration easier for patients and caregivers.



Ampule Accessories

[Order Online](#)[View More Products](#)

Ampule Opener

- Safer Design: Features a blade that scores the ampule stem to help ensure a clean break and reduce the chance of finger cuts.
- Compatible: Works with most ampules.
- Dimensions: 2½"L and 15/16" Dia.

Ampule Crackers

- Non-breakable plastic fracture and safety collar which is specifically designed to aid in the opening of pre-scored glass ampules
- Simply slip the Ampule Cracker over the top of the pre-scored portion, snap it, and discard both the ampule top and safety collar to prevent cross contamination
- Disposable: eliminates need for sterilization

Item No.	Description	Qty
7371	Ampule Opener	1
7372	Ampule Cracker for 1-4 ml Ampules	100
7372-10	Ampule Cracker for 1-4 ml Ampules	1000
7374	Ampule Crackers for 5-10 ml Ampules	100

TIP

Use an ampule opener or cracker to help achieve a clean, controlled break while reducing the risk of cuts and glass contamination.



Tube Set Fluid Transfer

Item No. 11

[Order Online](#)

[View More Products](#)

Features

- Sterile PVC tube set features silicone section
- Vented inlet spike and Luer lock outlet for filling sterile oral dispensers, Luer syringes, vials, elastomeric infusers, minibags and cassettes from IV bags or bottles of solution
- Not made with natural rubber latex
- **Quantity:** 10 per case



Admix Dispensing Pin

Item No. 350

[Order Online](#)

[View More Products](#)

Admix Dispensing Pins provide a simple, needle-free solution for the safe withdrawal and transfer of fluids from vials.

Features

- Sheathed plastic dispensing pin with 0.2 micron filtered vent (hydrophobic) and tethered cap.
- Use for vented additions and/or withdrawals of fluids from vials.
- Safe, simple, and needle-free
- **Dimensions:** 9.75 × 8.62 × 7.75"
- **Size:** 0.75" spike
- **Quantity:** 50 per package



TIP

Using a needle-free dispensing pin can help reduce the risk of accidental needle stick injuries and improve workplace safety.



Hypodermic Needles

[Order Online](#)[View More Products](#)

Designed for precision and reliability, hypodermic needles support comfortable injections, blood draws, and veterinary procedures.

Features

- Triple-bevel lubricated stainless steel bevel that minimizes patient discomfort
- Point durability optimized for vet procedure from vial to injection
- Ideal for blood draw
- Colour-coded hubs for easy needle gauge identification
- Packaged in soft packs for reduced waste
- Available in the most popular gauges and lengths



Item No.	Gauge	Needle Length	Qty
771810	18	1"	100
771815	18	1-1/2"	100
772010	20	1"	100
772015	20	1-1/2"	100
772034	20	3/4"	100
772234	22	3/4"	100
772210	22	1"	100
772506	25	5/8"	100

TIP

Color-coded needle hubs make it easy to identify the correct gauge at a glance, helping improve efficiency during veterinary procedures.



Needles & Syringe Combo

[Order Online](#)
[View More Products](#)

Designed for precision and reliability, hypodermic needles support comfortable injections, blood draws, and veterinary procedures.

Features

- Triple-bevel lubricated stainless steel bevel that minimizes patient discomfort
- Point durability optimized for vet procedure from vial to injection
- Ideal for blood draw
- Packaged in soft packs for reduced waste
- Available in the most popular gauges and lengths



Item #	Description	Volume	Needle Gauge	Needle Length
771001	Luer Slip	1mL	25G	5/8"
771010	Luer Slip	3mL	20G	1-1/2"
771009	Luer Slip	3mL	20G	1"
771013	Luer Slip	3mL	22G	1-1/2"
771012	Luer Slip	3mL	22G	1"
771011	Luer Slip	3mL	22G	3/4"
771002	Luer Lock	3mL	20G	1"
771003	Luer Lock	3mL	20G	1-1/2"
771004	Luer Lock	3mL	20G	3/4"
771005	Luer Lock	3mL	22G	3/4"
771006	Luer Lock	3mL	22G	1"
771007	Luer Lock	3mL	23G	1"
771008	Luer Lock	3mL	25G	5/8"

DID YOU KNOW?

Triple-bevel needles are designed to reduce penetration resistance, helping provide smoother injections and blood draws.



Adapter Caps

[Order Online](#)[View More Products](#)

Designed for use with liquid medication bottles and containers, adapter caps provide a secure, leak-resistant interface for oral dispensers, syringes, and other dispensing devices, helping ensure safe and accurate medication delivery.

Features

- Easy to Use: Converts a standard medication bottle into a filling device for plastic oral dispensers.
- Material: Low-density polyethylene

Item No.	Bottle Neck Size	Qty
8010-10	18 mm/A	20
8012-10	20 mm/B	20
7884-10	22 mm/C	20
7886-10	24 mm/D	20
8019-10	24 mm/M	20
7883-10	28 mm/E	20
8012-10	20 mm/B	20
8013-10	28 mm/F	20
8018-10	28 mm/K	20
8014-10	33 mm/G	20
19541	38 mm/I	20

TIP

Not sure which size to order?

Print the Adapter Cap Sizing Template to match the right adapter cap to your bottle. Print at 100% for the most accurate sizing.



Press-In Bottle Adapters

[Order Online](#)[View More Products](#)

Designed for use with liquid medication bottles, press-in bottle adapters create a secure, leak resistant connection for oral medication dispensers and syringes.

Features

- Mess-free filling: Easily fill oral dispensers without creating a mess
- Hassle-free design: Prescription bottle caps can be used with plugs in place. Adapters provide a no-leak closure
- Size options: Insert sizes correspond to the outer neck of bottle (not the inside opening)

Item No.	For Bottle Neck	Pack Qty	Case Qty
17282	20 mm	20	100
17283	24 mm	20	100
17284	28 mm	20	100
17285	33 mm	20	100
19392	22mm	20	100
19393	24mm	20	100
19394	28mm	20	100

Larger quantities available. Please contact us for more details.



TIP

Not sure which size to order?

Request a sample kit to find the best fit for your bottle.



SealSafe Slit Bottle Closures

[Order Online](#)[View More Products](#)


Features

- Prevent Contamination - Reseal quickly to provide a safe, clean and consistent seal. Helps prevent evaporation, spills and contamination.
- Easy Conversion - Non-sterile adapters quickly and easily convert standard prescription bottles into filling devices. Compatible with most brands of oral dispensers.
- Pierced Closure - Recommended for use with liquids slightly below, at or above the consistency of water.
- Cap Compatible - Original child-resistant and standard caps fit over closures.
- Material - Thermoplastic elastomer and a proprietary blend of polypropylene.

Item No.	Size	Pack Qty	Case Qty
20228	20 mm	20	200
20229	24 mm	20	200
20230	28 mm	20	200

TIP

Use bottle adapters to create a cleaner, more accurate liquid dispensing system while helping reduce spills, evaporation, and contamination.



Rapid Fill Connector

Luer Lock to Luer Lock with Cap

Item No. 13915

.....


[Order Online](#)

[View More Products](#)

Designed for accurate sterile filling applications, the Rapid Fill Connector creates a reliable bridge between luer connections and oral dispensers.

Features

- For filling sterile dispensers from a luer connection
- Luer lock-to-oral slip connection
- Connects to male luer or syringes, tube sets, repeating syringes, bag ports, etc.
- Flange facilitates contamination-free handling
- Sterile & latex free
- Made out of natural polypropylene
- Colour: Red
- Quantity: 50 per package, packaged individually



TIP

Use the Rapid Fill Connector's flange to maintain a secure grip and help reduce the risk of contamination during sterile transfer procedures.



PROTECT

Beard Guards

Masks


Gowns

Gloves

Infection Control

Spill Kits

Sharps Containers



Beard Covers & Bouffant Caps

Order Online

View More Products

Beard Covers

- Material - Polypropylene
- Sizes - Universal size
- Provides personal protection for doctors, surgeons, nurses and other medical professionals.

Beard Guards

- Completely cover beards, mustaches, and sideburns.
- Comfortable elastic strap loops behind ears

Tyvek® Bouffant Caps

- Material - Tyvek material offers superior breathability.
- Construction - Elastic headband provides an excellent barrier for dry particulates including microorganisms and more.



Item No.	Description	Qty
12077	Beard Covers	10
12077-31	Beard Covers	100
HB2060	Beard Guard	100
19346	Tyvek® Bouffant Caps	10
19347	Disposable SPP Bouffant Caps	100

DID YOU KNOW?

Proper hair and beard containment is an important part of contamination control, helping maintain a cleaner and more hygienic work environment.



Nitrile Medical Examination

[Order Online](#)

[View More Products](#)

Features:

- Material: Nitrile
- Thickness Indicator: Medium 5.0 mil
- Texture: Micro-textured
- Colour: Blue
- Available sizes: Small, Medium & Large. Specify when ordering

Item No.	Size	Qty
GNSMS	Small	100
GNSMM	Medium	100
GNSML	Large	100



Gowns

[Order Online](#)

[View More Products](#)

Disposable SMS Isolation Gown

- Disposable gown meets AAMI PB70 Level II barrier protection
- Gown made of non-woven material with full back closure.
- Features: Flat belt, elastic cuffs and tie-on neck closure.
- Colour: Yellow
- Size: Large

Item No.	Description	Qty
GW10Y	Disposable SMS Isolation Gown	10
GW100Y	Disposable SMS Isolation Gown	100



Face Masks

[Order Online](#)[View More Products](#)

3 Ply Face Mask with Ear Loop

- Elastic Ear Loops, latex free
- Disposable
- Pleated design for maximum comfort
- Lightweight & breathable
- Protection from non-toxic dust and contaminants
- Compliant to EN14683 Standard
- BFE≥98%, PFE ≥98%

Disposable Full Face Shield

- Construction - Foam headpiece with wide, soft, stretchable strap is lightweight and comfortable. Shield is anti-fogging.
- Dimensions - 7 1/2"L x 13"W

Ear Relief Bands

- Better Comfort - Relieve the tension and irritation caused by elastic bands that rub behind the ears after prolonged periods of face mask use
- Adjustable - Multiple notches on the semi-flexible plastic band allow users to choose their comfort level. Length can also be shortened using scissors.
- Dimensions - 8 5/16"L x 1 7/16"H



Item No.	Description	Qty
B344	3 Ply Face Mask with Ear Loop	50/box
12926	Disposable Full Face Shield	1
20445	Face Mask Ear Relief Bands	10

TIP

Amber sterile vials are Pair face masks with ear relief bands to reduce ear-loop strain and improve comfort during extended wear.



Spill Kits

[Order Online](#)
[View More Products](#)

Hazardous Drugs Spill Kit

- Rapid hazardous drug spill clean-up
- Enhanced staff and patient protection
- Convenient, ready-to-use emergency response solution
- Simplifies compliance with hazardous drug handling protocols

Emergency Spill Kit

- Complies with OSHA Blood borne Pathogens Standard.
- Includes
 - Stick-On Infectious Waste Disposal Bag, 12x14 (1)
 - Zip Closure Bag, 9x12 (1)
 - Pair of nitrile gloves (1)
 - P.A.W.S. Antimicrobial Hand Wipe (1)
 - High Absorbency Mop Pads (1)
 - Super Sani-Cloth Germicidal Disposable Wipe, XL (2)
 - Bilingual (English and Spanish) instruction sheet (1)

Kit Includes

- Nitrile Chemotherapy Gloves, Large* (2 pairs)
- Chemotherapy Gown* (1)
- Shoe Covers, non-slip (1 pair)
- Full-Face Safety Shield (1)
- 21g Green-Z® Fluid Solidifier (1)
- Scoop & Scraper (1)
- Yellow Chemotherapy Waste Bag, (1)
- Wiper Pads, 11" x 11" (2)
- Twist Tie (1)
- Identification Tag (1)
- Instruction Card* (1)
- Caution Sign (1)

Item No.	Description	Qty
20556	Hazardous Drugs Spill Kit	1
7464-01	Emergency Spill Kit	1

TIP

Place kits in high-risk areas and inspect contents regularly to ensure they are fully stocked and ready for use.



Chemo Plus Prep Mats

Item No. 9621-01

[Order Online](#)

[View More Products](#)

Features

- Material: The surface of the three-ply mat is firm and flat with a penetrable mesh top layer. The second layer is an absorbent gauze and the bottom layer is non-permeable plastic.
- Quantity: 50 per pack
- Dimensions: 17"W x 11"H



Sharps Containers

[Order Online](#)

[View More Products](#)

Features

- Material: Polypropylene
- Unused containers nest for easy storage
- Contains lock for final disposal
- Horizontal drop maximizes container space
- Safety mechanism prevents over-filling
- Prevents sharps from exiting the lid
- Bold red FULL indicator alerts user to replace container



Item No.	Capacity	Dimensions	Case Qty
9604-01	1 Quart	7.25" H x 6.5" W x 4.75" D	100
9601-01	5 Quart	10.75" W x 5.5" D	20

TIP

Replace the container when the red FULL indicator appears to help prevent overfilling and ensure safe sharps disposal.



Green-Z Solidifier, 15 oz

Order Online

View More Products

Overview

At times, workers may be exposed to blood and bodily fluid spills and other harmful substances such as toxic chemicals, biological hazards, and chemotherapy agents. To help reduce those exposure risks, rely on Green-Z® for improved safety in your spill cleanup operations.

Features:

- Green-Z® is chlorine-free and safe to use on chemotherapy, glutaraldehyde, and mercury spills.
- Contains fast-acting, super absorbent polymers.
- Chlorine-free (safe for use on chemotherapy, glutaraldehyde, and mercury spills).
- Eliminates splashing by solidifying and turning spills into semi-solid (gel) mass.
- Cost-effective cleanup solution. Dispose of infectious waste in accordance with federal regulations.
- Green-Z® is shelf-stable and does not require an expiration date. In order for the product to remain shelf-stable, it must be kept in a cool, dry place and the original packaging. If contents are hard or compact, simply shake the product to loosen before use.
- Packaging: Convenient shaker top bottle




TIP

Keep Green-Z® within easy reach to quickly solidify spills and minimize splashing, helping improve safety and simplify cleanup.



MONITOR

Temperature Thermometers
Cold Storage



Refrigeration/Freezer Thermometers

[Order Online](#)
[View More Products](#)

Features

- Triple display shows the current, minimum and maximum temperature reading, °C/°F switchable
- Unit displays exact time and date for all minimum and maximum temperature readings
- Glycol bottle probe insulates sensors from rapid temperature changes, ideal for mimicking vaccine temperatures
- Alarm alerts when temperature rises or falls outside set point
- Alarm is programmable in 0.1° increments; visual LED flashes and audible alarm sounds continuously, even if temperature returns to non-alarm range
- High-impact, chemical-resistant ABS plastic case
- For more information on each specific model, please visit the website [here](#)

Item No.	Description	Temp Range	Bottle Probe
10368	Traceable Memory Monitoring Thermometer	-40°C to 25°C	1-10ml glycol bottle
16088	Horizontal Thermometer	-50 to 70°C	-
17751	Traceable Jumbo Display Memory Monitoring Thermometer	-50 to 70°C	1-10ml glycol bottle
17752	Traceable Jumbo Display Memory Monitoring Thermometer	-50 to 70°C	1-5ml vaccine bottle
17988	Traceable Hi-Accuracy Refrigerator Thermometer	-50 to 70°C	1-10ml glycol bottle

DID YOU KNOW?

The included Traceable® calibration certificate is valid for one year. Many facilities replace their thermometers annually to maintain current compliance documentation.



Tenutō Fridges

[View All Cool Models](#)
[View Ultra Cold Models](#)
[View Combo Models](#)
[View Freeze Models](#)

Overview

Reliable cold storage solutions designed for healthcare, laboratory, research, veterinary, and scientific environments requiring consistent temperature performance and dependable storage protection.

From vaccines and medications to biological samples, reagents, and temperature-sensitive materials, professional environments require cold storage systems that deliver reliable and consistent performance every day. The Tenutō range of medical refrigerators and freezers is designed to support a wide range of applications across healthcare, laboratories, veterinary facilities, research institutions, pharmacies, and scientific environments. With solutions spanning cooling, freezing, combination, and ultra-low temperature storage, Tenutō medical coolers offers flexible configurations to support different operational and storage requirements.

Optimal for:

- Inpatient Pharmacy
- Outpatient Pharmacy
- Laboratory
- Nurse Ward
- Blood Depot



Item No.	Description	Temp Range	Capacity
21229	Cool 5L42 Countertop Medical Refrigerator	2°C~8°C	42 L / 1.48 Cu Ft
21300	Cool 5L126 Undercounter/Freestanding Medical Refrigerator	2°C~8°C	126 L / 4.45 Cu Ft
21301	Cool 5L316 Single Column Medical Refrigerator	2°C~8°C	316 L / 11.16 Cu Ft
22314	Cool 5L416 Single Column Medical Refrigerator	2°C~8°C	416 L / 14.69 Cu Ft

For more information, check out the specifications and benefits [here](#).

IDENTIFY & SECURE

Tamper-Evident Seals
Tape
Labels & Dispensers



Tamper Resistant Tape

[Order Online](#)
[View More Products](#)

Features

- **Convenient:** Perfect for use over items such as jar tops, bottle caps, syringes, boxes, specimen vials and similar containers that are hard to secure with ordinary security seals to help deter and show tampering. Tape molds easily to almost any item.
- **Deters Tampering:** Transparent Red polystyrene tape with serrated edges is printed with the message: "TAMPER RESISTANT TAPE" in Black and will tear in small pieces to show evidence of tamper upon removal.
- **Easy to Use:** Wax paper liner is slit lengthwise and allows for easy removal from roll.

Item No.	Description	Dimensions	Colour	Qty
12281	Integrity Tamper-Evident Tape Strips	6 $\frac{5}{8}$ " L x 1 $\frac{3}{8}$ " W	Red	100 Strips
8000	Tamper Resistant Tape	108' L x $\frac{1}{2}$ " H	Red	2 Rolls
8002	Tamper Resistant Tape	108' L x $\frac{1}{2}$ " H	Clear	2 Rolls
8005	Tamper Evident Tape, Trilingual	264" L x 2' H	Red	1 Roll
8007	Tamper Resistant Tape, French	108' L x 1" H	Red	1 Roll
8008	Tamper Tape	108' L x $\frac{1}{2}$ " H	Orange	2 Rolls
8009	Tamper Resistant Tape, French	108' L x $\frac{1}{2}$ " H	Orange	2 Rolls



DID YOU KNOW?

Flexible tamper resistant tape conforms easily to curved surfaces, making it ideal for bottles, jars, and other difficult-to-seal containers.



Tamper Evident Seals

[Order Online](#)[View More Products](#)

Tamper-Evident Seals

- Many Uses: Tamper-evident foil seals can also be used on all types of syringes, oral dispensers, vials, boxes, jars and more.
- Colors: Green, Light Blue, Red, Silver and Yellow.

Tamper-Evident Vial Seals

- Barrier Protection: Flexible foil seals offer a reliable moisture, particulate and bacterial barrier when properly applied.
- Sterile: Surface remains sterile until removed from the paper backing.
- Tamper-Evident: Needle punctures are highly visible in laminated foil.
- Compatible: 20mm bottle top
- Color Options: Blue, Green, Red and Yellow.



Item No.	Description	Dimensions	Qty
18812B	Tamper-Evident Seals	3-3/8"W x 1/2"H	1000 per roll
7867B	Tamper-Evident Vial Seals	20 mm	1000 per roll
7868	Tamper-Evident Syringe Seals	4 1/4"W x 7/16"H	1000 per roll

TIP

Apply tamper-evident seals to vials, syringes, and dispensers to help protect contents, provide visual evidence of tampering, and enhance medication safety.



Tape & Dispensers

[Order Online](#)
[View More Products](#)

Features


- Versatile: Transparent tape is designed for the protection and enhancement of prescription labels

Item No.	Tape Width	Tape Length	Core Size	Qty
T1-1	1" (24mm)	33 yards	1"	1 Roll
T1-3	1" (24mm)	72 yards	3"	1 Roll
T1-5	1.5" (38mm)	72 yards	3"	1 Roll
T1-7	2" (50mm)	55 yards	3"	1 Roll
TD1	Tape Dispenser For 1" Wide Tape (1" & 3" Cores)			1
GF2306	Tape Dispenser For 2" Wide Tape			1
GF2336	Tape Dispenser For 3" Wide Tape (3" Cores)			1







TIP


Individually packaged tape rolls help keep tape clean and ready for use until needed.



Administration/Instructions

			
LPH-58	2115	2284	8X
			
2X	6X	1-10G	1-354
			
1-436	1-44G	1-57G	P-1
			
4X	1-70G	15W	2285
			
1-60G	2287	1-669	1-49G
			
3W	10X	5X	2306
			
18620	2307	1-11G	1-389

Analgesics & Anti-Inflammatories

		
VET12	VET41	VET42





Anesthesia & Sedation

FENTANYL 2311	DEXMEDETOMIDINE _____ mg/mL CL-866	Methadone 10 mg/ml CL-1150	METHADONE LPH-123
Acepromazine _____ mg/mL VET26	Alfaxalone 10 mg/mL VET23	Diazepam 5 mg/mL VET11	Fentanyl 50 mcg/mL VET15
Ketamine 100 mg/mL VET24	Methadone _____ mg/mL VET36	Midazolam 5 mg/mL VET10	Propofol 10 mg/mL VET22
Remifentanyl _____ mg/mL VET16	Dexmedetomidine 500 mcg/mL VET27	Diazepam _____ mg/mL B2	Fentanyl _____ mcg/mL B25
Ketamine _____ mg/mL B41	PROPOFOL 17944	Remifentanyl _____ mcg/mL B46	VERSED®/Midazolam _____ mg/ _____ mL B39
Buprenorphine 0.5 mg/mL VET47	Buprenorphin 0.3 mg/mL VET35	Butorphanol 10 mg/mL VET40	Hydromorphone 2 mg/mL VET34
Hydromorphone 10 mg/mL VET33	Hydromorphone _____ mg/mL B12		

Antibiotics & Antimicrobials

CEFAZOLIN LPH-42	Enrofloxacin 50 mg/mL VET31	Metronidazole 5 mg/mL VET29	ANTIBIOTIC LPH-38
Ampicillin 100 mg/mL VET30	Cefazolin 100 mg/mL VET45		




Caution Labels

			
P-2	1-16G	AHSPT037	AHSPT07
			
2096	2805	8350	1-19G
			
18251	8353	FP-8353	18586
			
P-32			


Chemotherapy & Hazardous Drugs

			
18304	LPH-90	LPH-174	17520
			
LPH-56	2973	P-22	2373
			
PCT-1			

Communication Labels

			
P-11	2019	2215	2286
			
2156	1-13G	LPH-124	LPH-11
			
PP-2	11X	1-50G	1-38G
			
P-18			

Drug Identification Labels

			
LPH-193	LPH-46	LPH-136	CL-455
			
AHSPT014	VET25	B7	B27
			
B17	1-344	1-61G	

Emergency Drugs

Atropine 0.5 mg/mL	Glycopyrrolate 0.2 mg/mL	Atropine mg/mL	Glycopyrrolate mg/mL
VET13	VET14	B5	B6
norepinephrine	Dopamine 1600 mcg/mL	Epinephrine 1:1000	Norepinephrine mcg/mL
LPH-118	VET38	VET28	VET37
Ephedrine mg/mL	EPINEPHrine mg/mL	Naloxone mg/mL	Phenylephrine mg/mL
B31	B44	B28	B32

Expiration/Discard

1-35G	1-369	P-19	1-370
DATE OPENED: _____	Date Opened: _____ Date Expired: _____	DISCARD AFTER 28 DAYS EXP. DATE _____	EXPIRY DATE _____
P-21	2970	18369	P-26

GI & Supportive Care

ONDANSETRON	Famotidine 10 mg/mL	Metoclopramide 5 mg/mL	Ondansetron 2 mg/mL
LPH-34	VET32	VET43	VET49
Pantoprazole 4 mg/mL	Maropitant 10 mg/mL		
VET50	VET44		

Hormones & Endocrine

INSULIN	OXYTOCIN
2086	LPH-229

IV Solutions & Flushes

KCL ADDED	NORMAL SALINE	Hep-Saline	SALINE FLUSH
2304	LPH-64	LPH-21	2410
HEPARIN FLUSH	SALINE FLUSH		
17573	17574		

IV Solutions & Flushes

BUPIVACAINE 5 mg/mL	LIDOCAINE _____ mg/mL Date _____ Time _____ Int _____	Bupivacaine _____ mg/mL	Lidocaine 20 mg/mL
VET48	17793	LPH-87	VET19
Lidocaine mg/mL	ROPIVACAINE	Bupivacaine 5 mg/mL	Bupivacaine _____ mg/mL Date _____ Time _____ Int _____
B4	LPH-94	VET20	TP14

Refill Information

.....




CL-19



1-77G




2866



1-27G

Route of Administration

.....



17583



1-51G



1-78G



1-26G



1-18G



1-25G

Storage Information

.....



2083



1-15G



P-12



P-13



1-40G

Surgery/Recovery

.....

MEDICATION ADDED
DRUG ADDITIVE AMOUNT

ADDED BY _____ AM PM

DATE _____ TIME _____

PATIENT _____

ROOM NO. _____

START TIME _____

START DATE _____

FLOW RATE _____

BASE SOL'N _____

21-2

MEDICATION ADDED

DRUG ADDITIVE	AMOUNT

ADDED BY _____

DATE _____ TIME _____ A.M. P.M.

PATIENT _____

ROOM NO. _____

SPECIAL INFORMATION _____

MAR1

MEDICATION ADDED

PATIENT _____ RM. _____

DRUG _____

AMOUNT _____

ADDED BY _____

DATE _____ TIME _____

START TIME DATE . . . FLOW RATE

EXP. DATE _____

THIS LABEL MUST BE AFFIXED TO ALL INFUSION FLUIDS CONTAINING ADDITIONAL MEDICATION

MAR3

Write-on Labels, Blank

.....



2537



2536



2204



2553



2550



2211



A12



A10



A13



A11



2214



2552

Months of the Year, Package

.....



P-100
English



P-101
French

Label Dispensers & Accessories

[Order Online](#)

[View More Products](#)

Stack & Connect Label Dispenser

- **Connect Multiples:** Plastic label dispenser connects by the interlocking tabs on each side. Note: 10-roll dispensers do not connect to 20-roll dispensers.
- **#12721** includes nine removable dividers and accommodates roll and clamshell-packaged labels in quantities of 500 and 1,000.
- **#12998** does not include dividers. Accommodates roll and clamshell-packaged labels in quantities up to 500. Rolls of 1,000 labels will not work with this dispenser.
- **Holds Multiple Sizes:** Label slots will accommodate labels up to 5/8"W.
- **Mounting Options:** Mount under a counter or cabinet using Mounting Bracket (See options below).

Item No.	Description	Qty
12721	Stack & Connect Label Dispenser, 10-Roll	1
12998	Stack & Connect Label Dispenser, 20-Roll	1
17027	Mounting Bracket for #12721	1
17028	Mounting Bracket for #12998	1
18894	Peel N Go® Label Dispenser	1
19038	Mounting Bracket for Peel N Go® Label Dispenser	1
BK-10	Label Dispenser for 10 Rolls, Empty	1

TIP

Connect multiple dispensers and mount them under counters or cabinets to keep labels organized, accessible, and save valuable workspace.



ORDER FORM

Item No.	Description	Quantity

Pharmacy Name:		Contact Name:	
Address:		Phone:	
City:		Fax:	
Province:	Postal Code:	Email:	
Date:	Wholesaler:	Account No.:	

.....

Please email or fax the completed order form to us.

orders@pharmasystems.com

fax: 1.888.475.7155

Please note: Shipping charges will be added to each order

 **PharmaSystems®**

pharmasystems.com