Introducing

Tenutō Cool Pharmacy Refrigerators

Purpose-Built Refrigerators that Help Protect High-Value Medications and Critical Inventory





Tenutō Cool refrigerators ensure precise temperature control, unmatched reliability, and seamless workflow integration. From LBA foaming technology to electric heating and e-glass, Tenutō Cool has all the features your facility needs for trusted cold chain storage.

Stable Temperature Control
Designed to maintain reliable cooling performance for temperature-sensitive materials stored between 2°C to 8°C.

Practical Interior Design
Multiple storage capacities, and adjustable shelving & organized storage layouts help support accessibility and efficient space utilization.

Monitoring & Alarm Functions
Select models include temperature log, monitoring, visual and audible high and low temperature alarms, and operational safety features to support storage protection.

Professional Design for Critical Environments
Designed for use across healthcare applications requiring dependable refrigerated storage and cold chain management.

Ready for Every Stage of Your Pharmacy's Growth

- Replace Aging Units
- Expand Cold Storage Capacity
- Equip New Locations Opening

Now Available in Canada, Exclusively from PharmaSystems

Trusted Partner
Your single source for pharmacy supplies. Servicing Canadian pharmacies since 1977.

Easy Ordering
Quick & easy ordering online, or by email, phone or fax, with same-day order processing.

Dedicated Support
Connect with our bilingual team of experts for product support and recommendations.

PharmaSystems®

PharmaSystems.com/Tenuto
PharmaSystems Inc. 151 Telson Road, Markham, Ontario
Email: info@PharmaSystems.com | **Call:** 1.888.475.2500

Tenutō Cool Models

Countertop, Under-counter, and Single-Column Models Designed for Reliable Cold Storage Performance

Click each model to view specifications, pricing, and order online.







		21299	21300	21301
Performance	**Capacity**	42 L / 1.48 Cu Ft	126 L / 4.45 Cu Ft	316 L / 11.16 Cu Ft
	Type	Countertop	Undercounter	Single Column
	Interior Dimension (W*D*H)	14.96×11.81×16.14 in	19.88×18.39×24.41 in	19.69×18.27×52.87 in
	Exterior Dimension (W*D*H)	17.72×18.07×22.42 in	23.43×24.21×31.89 in	23.43×23.74×75.59 in
	Top Cooling Technology	N	N	Y
	Temperature Range (°C)	2°C~8°C	2°C~8°C	2°C~8°C
	Accuracy (°C)	0.1	0.1	0.1
	Fluctuation(°C)	≤3	≤2	≤2
	Uniformity(°C)	±1	±1	±1
Functions	**High/Low Temperature Alarm**	✅	✅	✅
	Power Failure Alarm	✅	✅	✅
	Sensor Error Alarm	✅	✅	✅
	Low Battery Alarm	✅	-	✅
	High Ambient Temp Alarm	✅	✅	✅
	Door Ajar Alarm	✅	✅	✅
	Communication failure Alarm	✅	-	✅
	Recorder fault alarm	✅	-	✅
	Recording is not started Alarm	✅	-	✅
	Max and min temperature check	✅	✅	✅
	Dual Display (temperature & humidity)	-	-	✅
	Inner Light	LED Light	LED Light	LED Light
Accessories	**Caster**	-	-	4
	Leveling Feet	4	4	-
	Door Lock	Standard	Standard	Standard
	Shelves/Drawers	1/1	3/-	5/-
	Porthole (Qty/Diameter mm)	1/25mm (0.98 in)	1/25mm (0.98 in)	1/25mm (0.98 in)
	USB Interface	Standard	Standard	Standard
	Remote Alarm Contact	Standard	Standard	Standard
	Power Failure Backup System (Alarm)	16h	20h	16h
Technical Data	**Cooling Type/Defrost Mode**	Forced Air Cooling/Auto	Forced Air Cooling/Auto	Forced Air Cooling/Auto
	Compressor	GMCC	GMCC	GMCC
	Refrigerant Gas	R600a	R600a	R600a
	Noise Level (dB(A))	45	45	220-240V/50Hz: 52 115V/60Hz: 60 220-240V/60Hz: 52
	Certifications	CE, UL, NSF, Energy Star	CE, UL, NSF, Energy Star	CE, UL, NSF, Energy Star
	Controller	Microprocessor	Microprocessor	Microprocessor
	Display	LED	LED	LED

PharmaSystems®

PharmaSystems.com/Tenuto
PharmaSystems Inc. 151 Telson Road, Markham, Ontario
Email: info@PharmaSystems.com | **Call:** 1.888.475.2500